
PreON User Manual

Document Information

SUBJECT	DESCRIPTION
PRODUCT:	PreON
DOCUMENT TYPE:	User Manual
DOCUMENT NAME:	PreON User Manual
DOCUMENT NUMBER:	UM-PN-1001
REVISION:	06
DATE:	06/05/2022
COMPANY:	PreOmics GmbH

Document Change Control

Adaptations and updates of the present document are a consequence of technical improvements and amendments of the PreON. Any changes within this document are to be controlled.

Document revisions

REVISION	DOCUMENT FILE NAME	REASON FOR CHANGE	BY	VALID FROM
01	SM-PN1001-01	New Document	HSE	08/10/2019
02	SM-PN1001-02	Updates for SW 1.2	HSE	20/02/2020
03	SM-PN1001-03	Updates for SW 1.4	HSE	13/10/2020
04	SM-PN1001-04	Updates for SW 1.4	HSE	29/01/2021
05	SM-PN1001-05	Updates for PreON fuses	HSE	22/03/2021
06	SM-PN1001-06	Updates for SW 1.10	HSE	06/05/2022

Trademarks

Trademarks: PreON

© 2020 PreOmics GmbH, all rights reserved.

Table of content

1	Introduction	8
1.1	About this user manual.....	8
1.1.1	Technical assistance change	9
1.1.2	Policy statement.....	9
1.2	Intended use of the PreON.....	9
1.3	Requirements for PreON users	9
2	Safety Information	10
2.1	Proper use	10
2.2	Electrical safety.....	12
2.3	Environment.....	13
2.3.1	Operating conditions	13
2.4	Biological safety.....	13
2.4.1	Samples.....	13
2.5	Chemicals	14
2.5.1	Toxic fumes	14
2.6	Waste disposal.....	14
2.7	Mechanical hazards.....	15
2.7.1	Centrifuge.....	15
2.8	Heat hazard	16
2.9	Maintenance safety.....	16
2.10	Radiation safety	17
2.11	Symbols on the PreON.....	18
3	General Description	19
3.1	External features of the PreON.....	21
3.2	Internal features of the PreON	25
4	Installation Procedures.....	32
4.1	Installation environment.....	32
4.1.1	Site requirements	32
4.1.2	Power requirements.....	33
4.1.3	Grounding requirements	33
4.2	Installation of AC power cord.....	34
4.3	Unpacking the PreON.....	34
4.4	Installing the PreON.....	35
4.4.1	Remove PreON accessories.....	35
4.4.2	Remove the foam protector.....	35

4.4.3	Installation of AC power cord.....	35
4.4.4	Powering ON the PreON	36
4.4.5	Installing the centrifuge rotor and buckets	36
4.5	Configuration of the PreON	37
4.5.1	Setting basic system data	38
4.5.2	User settings	39
4.5.3	Managing users.....	40
4.5.4	Changing password	43
4.5.5	Managing protocols.....	44
4.5.6	Updating software.....	46
4.6	Packing the PreON.....	48
5	Operating Procedures	50
5.1	PreON software	53
5.1.1	Entering text and numbers.....	54
5.2	Switching on the instrument and logging in	55
5.3	Setting up a protocol run.....	56
	iST-NHS protocol.....	58
	Step 1 Scan barcode chemical labeling reagent	58
	Step 2 Select sample type	58
	Step 3 Select digestion duration.....	58
	Step 4 Select number of samples	59
	Step 5 Activate Low Input protocol	59
	Step 6 Select chemical labeling plex degree	59
	Step 7 Select amount of chemical labeling reagent	60
	Step 8 Chemical labeling options.....	60
	Step 9 Select labeling time	60
	Step 10 Select split or complete protocol	61
	Step 11 Load ELUTE later.....	61
	Step 12 Sample drying	61
	Step 13 Review Settings	62
	Step 14 Insert 1.....	62
	Step 15 Insert 3.....	62
	Step 16 Load reagents	63
	Step 17 Load tip racks.....	63
	Step 18-1 Load the centrifuge rotor adapter – Spin Column.....	63
	Step 18-2 Load the centrifuge rotor adapter – Middle Position.....	64

Step 18-3 Load the centrifuge rotor adapter – Elute Position.....	64
Step 18-4 Load the centrifuge rotor adapter.....	64
Step 19 Load shaker.....	65
iST protocol	66
Step 1 Select sample type	66
Step 2 Select digestion duration.....	66
Step 3 Select number of samples	66
Step 4 Low input protocol.....	67
Step 5 Load ELUTE later.....	67
Step 6 Sample drying.....	67
Step 7 Review	68
Step 8 Insert 1	68
Step 9 Insert 3	68
Step 10 Load reagents	69
Step 11 Load tip racks.....	69
Step 12-1 Load the centrifuge rotor adapter – Spin Column.....	69
Step 12-2 Load the centrifuge rotor adapter – Middle Position.....	70
Step 12-3 Load the centrifuge rotor adapter – Elute Position.....	70
Step 12-4 Load the centrifuge rotor adapter.....	70
Step 13 Load shaker.....	71
iST-BCT protocol	72
Step 1 Select sample type	72
Step 2 Select digestion duration.....	72
Step 3 Select number of samples	72
Step 4 Low input protocol.....	73
Step 5 Load ELUTE later.....	73
Step 6 Sample drying.....	73
Step 7 Review	74
Step 8 Insert 1	74
Step 9 Insert 3	74
Step 10 Load reagents	75
Step 11 Load tip racks.....	75
Step 12-1 Load the centrifuge rotor adapter – Spin Column.....	75
Step 12-2 Load the centrifuge rotor adapter – Middle Position.....	76
Step 12-3 Load the centrifuge rotor adapter – Elute Position.....	76
Step 12-4 Load the centrifuge rotor adapter.....	76

Step 13 Load shaker	77
5.3.1 Loading reagents and tips	77
5.3.2 Loading the centrifuge.....	78
5.3.3 Loading the shaker	81
5.4 Starting the protocol run	82
5.5 End of the protocol run	83
5.6 Stopping a protocol run.....	84
5.7 Saving run reports to the USB stick.....	84
5.8 Logging out.....	85
6 Maintenance Procedures	85
6.1 Cleaning agents	86
6.1.1 Disinfection	86
6.2 Servicing	87
6.3 Regular maintenance procedure	88
6.4 Daily maintenance procedure.....	89
6.5 Monthly maintenance procedure.....	90
6.6 Periodic maintenance procedure	90
6.6.1 Cleaning the robotic arm modules	91
6.6.2 Cleaning the centrifuge	91
6.6.3 Tightness test.....	95
6.6.4 O-Ring exchange	96
6.7 Decontaminating the PreON	100
7 Troubleshooting	100
7.1 General information.....	100
7.2 Contacting PreOmics Technical Services.....	100
7.2.1 Creating a support package	101
7.3 Operation	101
7.3.1 Protocol interruption	103
7.3.2 Centrifuge.....	103
7.3.3 Reagent volume detection and ultrasonic pipe.....	105
7.3.4 Touchscreen.....	105
8 Automated Status Supervision.....	108
8.1 HTML Website	108
8.2 /status.ini file	108
8.2.1 Format.....	109
9 Glossary.....	110

10	Appendix A – Technical data.....	112
10.1	Operating conditions	112
10.2	Transport conditions.....	112
10.3	Storage conditions	112
10.4	Mechanical data and hardware features	112
10.5	Declaration of Conformity	114
10.6	Waste Electrical and Electronic Equipment (WEEE)	114
10.7	FCC Declaration	114
10.8	Liability Clause.....	115
11	Appendix B – PreON Accessories	116
12	Appendix C – Consignes de sécurité.....	117
12.1	Utilisation appropriée	117
12.2	Sécurité électrique	119
12.3	Atmosphère.....	120
12.4	Sécurité biologique.....	120
12.4.1	Échantillons	120
12.5	Produits chimiques.....	121
12.5.1	Vapeurs toxiques.....	121
12.6	Traitement des déchets.....	122
12.7	Dangers mécaniques.....	122
12.7.1	Centrifuge.....	122
12.8	Dangers liés à la chaleur	123
12.9	Sécurité relative à la maintenance	124
12.10	Sécurité contre les rayonnements	125
12.11	Symboles sur le PreON.....	125
13	Appendix D – Sicherheitshinweise	126
13.1	Sachgemäße Handhabung	126
13.2	Schutz vor Stromschlag	128
13.3	Atmosphären.....	129
13.3.1	Operating conditions	129
13.4	Biologische Sicherheit.....	130
13.4.1	Proben.....	130
13.5	Chemikalien.....	131
13.5.1	Giftige Dämpfe	131
13.6	Entsorgen von Abfällen	131
13.7	Gefahren durch mechanische Teile	132

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1 Introduction

Thank you for choosing the PreON. We are confident it will become an integral part of your laboratory. Before using the PreON, it is essential that you read this user manual carefully and pay attention to the safety information. The instructions and safety information in the user manual must be followed to ensure safe operation of the instrument and to maintain the instrument in a safe condition.

Only professionals who understand the user manual in the introduction and safety information section must operate the device.

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1.1 About this user manual

This user manual provides information about the PreON in the following sections:

1. Introduction
2. Safety Information
3. General Description
4. Installation Procedures
5. Operating Procedures
6. Maintenance Procedures
7. Troubleshooting
8. Automated Status Supervision
9. Appendix A – Technical data
10. Appendix B – PreON Accessories
11. Appendix C – Consignes de sécurité (Safety Information in French)
12. Appendix D – Sicherheitshinweise (Safety Information in German)
13. Appendix E – Informazioni di sicurezza (Safety Information in Italian)

1.1.1 Technical assistance change

At PreOmics, we pride ourselves on the quality and availability of our technical support. Our Technical Services Departments are staffed by experienced scientists with extensive practical and theoretical expertise in molecular biology and the use of PreOmics products. If you have any questions or experience any difficulties regarding the PreON or PreOmics products in general, do not hesitate to contact us.

PreOmics customers are a major source of information regarding advanced or specialized uses of our products. This information is helpful to other scientists as well as to the researchers at PreOmics. We therefore encourage you to contact us if you have any suggestions about product performance or new applications and techniques.

For technical assistance and more information, please call one of the PreOmics Technical Service Departments or local distributors (see back cover or visit www.PreOmics.com).

1.1.2 Policy statement

It is the policy of PreOmics to improve products as new techniques and components become available. PreOmics reserves the right to change specifications at any time.

To produce useful and appropriate documentation, we appreciate your comments on this user manual. Please contact PreOmics Technical Services.

1.2 Intended use of the PreON

The PreON is designed to perform fully automated purification of peptides and proteins in molecular biology applications. The system is intended for use by professional users, such as technicians and physicians trained in molecular biological techniques and the operation of the Proton. The Proton system is not intended for diagnostic use.

The PreON is intended to be used only in combination with PreOmics kits indicated for use with the PreON for the applications described in the kit handbooks.

1.3 Requirements for PreON users

The table below covers the general level of competence and training necessary for transportation, installation, use, maintenance, and servicing of the PreON.

Task	Personnel	Training and experience
Transportation	No special requirements	No special requirements
Storage	No special requirements	No special requirements
Routine use and maintenance	Laboratory technicians or equivalent	Professional personnel, such as technicians and physicians, who are trained in molecular biological techniques and in the operation of the proton.
Installation and servicing	PreOmics Field Service Specialists only	Trained and authorized by PreOmics

2 Safety Information

Before using the PreON, it is essential that you read this user manual carefully and pay attention to the safety information. The instructions and safety information in the user manual must be followed to ensure safe operation of the instrument and to maintain the instrument in a safe condition. **Note:** Translations of the Safety Information in French, German and Italian are available in Appendix C – Consignes de sécurité, Appendix D – Sicherheitshinweise and Appendix E – Informazioni di sicurezza.

The following types of safety information appear in this manual.

WARNING

The term WARNING is used to inform you about situations that could result in personal injury to you or other persons.
Details about these circumstances are given in a box like this one.

CAUTION

The term CAUTION is used to inform you about situations that could result in damage to the instrument or other equipment.
Details about these circumstances are given in a box like this one.

Important note: All symbols that are marked on the PreON are explained in detail in the related sections of this user manual. In particular it is mandatory to consult the user manual in case of the generic symbol shown below and all specific symbols that are safety related.

Safety related information. User manual must be consulted to find out the nature of potential hazards and any actions which have to be taken to avoid them

The advice given in this manual is intended to supplement, not supersede, the normal safety requirements prevailing in the user's country.

2.1 Proper use

WARNING

Risk of personal injury and material damage

[W1]

Improper use of the PreON may cause personal injuries or damage to the instrument. The PreON must only be operated by professional personnel who have been appropriately trained. Servicing of the PreON must only be performed by a PreOmics Field Service specialist.

Perform the maintenance as described in section Maintenance Procedures. PreOmics charges for repairs that are required due to incorrect maintenance.

WARNING

Risk of personal injury and material damage

[W2]

The PreON is too heavy to be lifted by one person. To avoid personal injury or damage to the instrument, do not lift the instrument alone.

WARNING

Risk of personal injury and material damage

[W3]

Do not attempt to move the PreON during operation.

CAUTION

Damage to the instrument

[C1]

Avoid spilling water or chemicals onto the PreON. Damage caused by water or chemical spillage will void your warranty.

In case of emergency, power OFF the PreON at the power switch located in front of the instrument and unplug the power cord from the power outlet.

CAUTION

Damage to the instrument

[C2]

Only use PreOmics spin column-based cartridges and PreON specific consumables with the PreON. Damage caused by use of other types of spin column-based cartridges or chemistries will void your warranty.

WARNING

Risk of personal injury and material damage

[W4]

Do not use damaged rotor adapters. The rotor adapters can only be used once. High g forces exerted in the centrifuge can cause damage to re-used rotor adapters.

CAUTION

Damage to the instrument

[C3]

Empty the tip disposal container prior to use to prevent a tip jam in the waste drawer. Failure to empty the waste container may block the robotic arm that could cause run failure or instrument damage.

WARNING

Risk of personal injury and material damage

[W5]

To prevent plastic crash, load the tubes properly. After a plastic crash, sharp plastic particles could be inside the centrifuge. Be careful when handling items inside the centrifuge.

CAUTION

Damage to the instrument

[C4]

Only use the correct volume of liquids.

Exceeding the recommended volume of liquids may damage the centrifuge rotor or instrument.

WARNING

Risk of fire or explosion

[W6]

When using flammable liquids on the PreON (Acetonitrile, Isopropyl alcohol), handle such liquids carefully and in accordance with the required safety regulations. If liquid has been spilled, wipe it off and leave the PreON hood open to allow flammable vapors to disperse.

WARNING

Risk of fire or explosion

[W7]

The PreON is intended for use with reagents and substances supplied with PreOmics kits. Use of other reagents and substances may lead to fire or explosion.

If hazardous material is spilled on or inside the PreON, the user is responsible for carrying out appropriate decontamination.

CAUTION

Damage to the instrument

[C5]

Avoid contact of the reagents with the instrument surfaces and centrifuge (rotor, buckets, gasket) as this will damage instrument surfaces and centrifuge gasket and cause corrosion of metal on the worktable, the centrifuge rotor and buckets.

Clean spillage on the worktable and centrifuge immediately after the run.

Note: Do not place items on top of the PreON hoods.

CAUTION

Damage to the instrument

[C6]

Do not lean against the touchscreen when it is pulled out.

2.2 Electrical safety

Note: Disconnect the line power cord from the power outlet before servicing.

WARNING

Electrical hazard

[W8]

Any interruption of the protective conductor (earth/ground lead) inside or outside the instrument or disconnection of the protective conductor terminal is likely to make the instrument dangerous.

Intentional interruption is prohibited.

Lethal voltages inside the instrument

When the instrument is connected to line power, terminals may be live and opening covers or removing parts is likely to expose live parts.

WARNING

Damage to electronics

[W9]

Before powering ON the instrument, make sure that the correct supply voltage is used.

Use of incorrect supply voltage may damage the electronics.

To check the recommended supply voltage, refer to the specifications indicated in the type plate of the instrument.

WARNING

Risk of electric shock

[W10]

Do not open any panels on the PreON.

Risk of personal injury and material damage may occur.

Only perform maintenance that is specifically described in this user manual.

To ensure satisfactory and safe operation of the PreON, follow these guidelines:

- The line power cord must be connected to a line power outlet that has a protective conductor (earth/ground).
- Do not adjust or replace internal parts of the instrument.
- Do not operate the instrument with any covers or parts removed.

If the instrument becomes electrically unsafe, prevent other personnel from operating it, power OFF the instrument, disconnect it from the power outlet and contact PreOmics Technical Services.

The instrument may be electrically unsafe when:

- It or the line power cord appears to be damaged.
- It has been stored under unfavorable conditions (see chapter 9.3 Storage conditions) for a prolonged period.
- It has been subjected to severe transport stresses. Indications can be that the transportation crate is damaged, or the instrument housing is compromised.
- Liquids come in contact directly with electrical components of the PreON.

2.3 Environment

2.3.1 Operating conditions

WARNING

Explosive atmosphere

[W11]

The PreON is not designed for use in an explosive atmosphere.

CAUTION

Damage to the instrument

[C7]

Direct sunlight may bleach parts of the instrument and cause damage to plastic parts.

The PreON must be located out of direct sunlight.

2.4 Biological safety

Specimens and reagents containing materials from humans should be treated as potentially infectious. Use safe laboratory procedures as outlined in publications such as Biosafety in Microbiological and Biomedical Laboratories, HHS (www.cdc.gov/labs/pdf/CDC-BiosafetyMicrobiologicalBiomedicalLaboratories-2009-P.PDF).

2.4.1 Samples

Samples may contain infectious agents. You should be aware of the health hazard presented by such agents and should use, store, and dispose of such samples according to the required safety regulations.

WARNING

Samples containing infectious agents

[W12]

Some samples used with this instrument may contain infectious agents. Handle such samples with the greatest of care and in accordance with the required safety regulations.

Always wear safety glasses, 2 pairs of gloves, and a lab coat.

The responsible body (e.g., laboratory manager) must take the necessary precautions to ensure that the surrounding workplace is safe, and that the instrument operators are suitably trained and not exposed to hazardous levels of infectious agents as defined in the applicable Safety Data Sheets (SDSs) or OSHA¹ or ACGIH² or documents.

Venting for fumes and disposal of wastes must be in accordance with all national, state, and local health and safety regulations and laws.

¹ OSHA : Occupational Safety and Health Administration (United States of America).

² ACGIH : American Conference of Government Industrial Hygienists (United States of America).

2.5 Chemicals

WARNING

Hazardous chemicals

[W13]

Some chemicals used with this instrument may be hazardous or may become hazardous after completion of the protocol run.

Always wear safety glasses, gloves, and a lab coat.

The responsible body (e.g., laboratory manager) must take the necessary precautions to ensure that the surrounding workplace is safe and that the instrument operators are not exposed to hazardous levels of toxic substances (chemical or biological) as defined in the applicable Safety Data Sheets (SDSs) or OSHA or ACGIH documents.

Venting for fumes and disposal of wastes must be in accordance with all national, state, and local health and safety regulations and laws.

CAUTION

Damage to the instrument

[C8]

Avoid contact of the reagents with the instrument surfaces and centrifuge (rotor, buckets, gasket) as this will damage instrument surfaces and centrifuge gasket and cause corrosion of metal on the worktable, the centrifuge rotor and buckets.

Clean spillage on the worktable and centrifuge immediately after the run.

2.5.1 Toxic fumes

If working with volatile solvents or toxic substances, you must provide an efficient laboratory ventilation system to remove vapors that may be produced.

WARNING

Vapors from volatile solvents

[W14]

Some of the solvents used with the instrument are volatile. Use the instrument only with adequate ventilation according to the applicable Safety Data Sheets (SDSs) or OSHA or ACGIH[†] documents.

WARNING

Toxic fumes

[W15]

Do not use bleach to clean or disinfect the PreON. Bleach in contact with salts from the buffers can produce toxic fumes.

WARNING

Toxic fumes

[W16]

Do not use bleach to disinfect used labware. Bleach in contact with salts from the buffers used can produce toxic fumes.

2.6 Waste disposal

Used labware, such as sample tubes, PreOmics spin column-based cartridges, filter-tips, buffer bottle and enzyme tubes, or rotor adapters, may contain hazardous chemicals or infectious agents from the purification process. These hazardous wastes must be collected and disposed of properly according to local safety regulations.

For more information about how to dispose of the PreON, see section 10.6 Waste Electrical and Electronic Equipment (WEEE).

WARNING

Hazardous chemicals and infectious agents

[W17]

The waste may contain toxic material and must be disposed of properly. Refer to your local safety regulations for proper disposal procedures.

2.7 Mechanical hazards

The hood of the PreON must remain closed during operation of the instrument. Only open the hood when instructed to do so by the instruction for use.

While loading the worktable, always stand clear of the instrument. Do not lean on the worktable when the robotic arm of the instrument is moving to reach loading position with its lid open. Wait until the robotic arm completed its movements before you start to load or unload.

WARNING

Moving parts

[W18]

Avoid contact with moving parts during operation of the PreON. Under no circumstances should you place your hands under the robotic arm when it is lowering. Do not attempt to move any tip racks or tubes whilst the instrument is operating.

WARNING

Moving parts

[W19]

To avoid contact with moving parts during operation of the PreON, the instrument must be operated with the hood closed.

If the hood sensor or lock is not functioning correctly, contact PreOmics Technical Services.

2.7.1 Centrifuge

Make sure that the rotor and buckets are installed correctly. All buckets must be mounted before starting a protocol run, regardless of the number of samples to be processed. If the rotor or buckets show signs of mechanical damage or corrosion, do not use the PreON; contact PreOmics Technical Services.

CAUTION

Damage to the instrument

[C9]

The PreON must not be used if the centrifuge lid is broken, or if the lid lock is damaged.

Make sure that no loose material is inside the centrifuge during operation.

Make sure that the rotor is installed correctly and that all buckets are properly mounted, regardless of the number of samples to be processed. Load the rotor only as instructed by the software.

Only use rotors, buckets, and consumables designed for use with the PreON. Damage caused by use of other consumables will void your warranty.

We recommend replacing the centrifuge after 20,000 cycles, which is equivalent to 9 years of usage with two runs per day for 220 days each year. For more information contact PreOmics Technical Services.

In case of breakdown caused by power failure, the centrifuge lid can be opened manually to remove the samples (see section 7.3.2).

Cleaning and decontamination may be necessary as a safeguard before LABORATORY, CENTRIFUGES, ROTORS, and any accessories are maintained, repaired, or transferred. Manufacturers may provide a format for users to document that such treatment has been carried out. For more information contact PreOmics Technical Services.

WARNING

Moving parts

[W20]

In case of breakdown caused by power failure, remove the power cord and wait 10 minutes before attempting to manually open the centrifuge lid.

CAUTION

Damage to the instrument

[C10]

After a power failure, do not move the z-module (robotic arm) manually. Damage may occur if the centrifuge lid is closed and collides with the z-module.

WARNING

Risk of personal injury and material damage

[W21]

Raise the centrifuge lid carefully. The lid is heavy and may cause injury if it falls.

CAUTION

Risk of overheating

[C11]

To ensure proper ventilation, maintain a minimum clearance of 10 cm at the sides and rear of the PreON.

Slits and openings that ensure the ventilation of the PreON must not be covered.

2.8 Heat hazard

The PreON worktable contains a heated shaker.

WARNING

Hot surface

[W22]

The shaker can reach temperatures of up to 83°C (181.4°F). Avoid touching it when it is hot. Carefully remove the samples after a run.

2.9 Maintenance safety

**WARNING/
CAUTION**

Risk of personal injury and material damage

[W23]

Only perform maintenance that is specifically described in this user manual.

WARNING

Risk of explosion

[W24]

When cleaning the PreON with alcohol-based disinfectant, leave the PreON hood open to allow flammable vapors to disperse.

Only clean the PreON when worktable components have cooled down.

WARNING 	Risk of fire [W25] Do not allow cleaning fluid or decontamination agents to come into contact with the electrical parts of the PreON.
WARNING 	Risk of personal injury and material damage [W26] To prevent the rotor nuts from loosening during operation of the centrifuge, securely tightened the nuts using the rotor key supplied with PreON.
WARNING 	Risk of personal injury and material damage [W27] Make sure that lids from spin column-based cartridges and 1.5 ml microcentrifuge tubes are in the correct position and pushed all the way down to the bottom of the slots on the sides of the rotor adapter. Incorrectly positioned lids can break off during centrifugation.
WARNING 	Risk of personal injury and material damage [W28] Be sure the lid is completely removed from the spin column. Spin column-based cartridges with partially removed lids may not be removed properly from the rotor, causing the protocol run to crash.
CAUTION 	Damage to the instrument [C12] Do not use bleach, solvents, or reagents containing acids, alkalis or abrasives to clean PreON.
CAUTION 	Damage to the instrument [C13] Avoid contact of the reagents with the instrument surfaces and centrifuge (rotor, buckets, gasket) as this will damage instrument surfaces and centrifuge gasket and cause corrosion of metal on the worktable, the centrifuge rotor and buckets. Clean spillage on the worktable and centrifuge immediately after the run.
CAUTION 	Damage to the instrument [C14] Do not use spray bottles containing alcohol or disinfectant to clean surfaces of the PreON. Spray bottles should be used only to clean items that have been removed from the worktables.

2.10 Radiation safety

WARNING 	Risk of personal injury [W29] Hazard Level 2 laser light: Do not stare into the light beam when using handheld bar code scanner, Read for further instructions the user manual of the handheld scanner.
---	---

2.11 Symbols on the PreON

Symbol	Location	Description
	Next to the shaker	Hot surface – the temperature of the shaker can reach up to 83°C (181.4°F) (cf. section 2.8)
	Near the centrifuge; near the robotic arm	Safety related information. User manual must be consulted to find out the nature of potential hazards and any actions which have to be taken to avoid them Mechanical hazard – avoid contact with moving parts (cf. section 2.7).
	On the instrument, near bottle rack	Safety related information. User manual must be consulted to find out the nature of potential hazards and any actions which have to be taken to avoid them Fire hazard – usage of ethanol in bottle rack (cf. section 2.1).
	In front of worktable Inside the waste drawer	Biohazard risks – substances in the waste drawer or used with this instrument may be hazardous. based on biological activity must be handled with gloves (cf. section 2.4).
	Back of the instrument next to the mains connection	Possibility of electric shock. Do not open any panels on the PreON (cf. section 2.2)
	Front of Barcode- scanner	Laser radiation - do not stare into the light beam (cf. 2.10)
	Type plate on the back of the instrument	CE marking to indicate conformity with the applicable requirements of the European Community
	Type plate at the back of the instrument	CSA marking to indicate compliance to applicable CSA and ANSI/UL Standards, for use in Canada and the U.S..
	Type plate on the back of the instrument	Waste Electrical and Electronic Equipment (WEEE) marking for Europe to indicate that product must not be disposed of with other waste (cf. section 9.6)
	Type plate on the back of the instrument	Manufacturer
	Type plate on the back of the instrument	Serial number of the instrument

3 General Description

The PreON performs fully automated processing of up to 12 samples (up to 16 samples with iST-NHS and TMT16plex) per run. The PreON is designed to automate selected PreOmics kits. The PreON controls integrated components, including a centrifuge, heated shaker, pipetting system and robotic gripper.

The PreON is preinstalled with iST, iST-NHS and iST-BCT protocols for processing PreOmics spin column-based cartridges for sample preparation and peptide cleanup. The user selects a protocol using the touchscreen and loads labware, samples, and reagents onto the PreON worktable. The user then closes the instrument hood and starts the protocol, which provides all necessary commands for sample lysis and purification using PreOmics spin column-based cartridges. A fully automated load check helps to ensure correct loading of the worktable.


The range of protocols available is continually expanding.

PreON principle


Sample preparation using the PreON follows the same steps as the manual procedure (i.e., lyse, proteolytically digest, wash and elute). No change in purification chemistry is required as you simply continue to use trusted PreOmics kits.

1. Samples are lysed in the heated orbital shaker.
2. Each lysate is digested with Proteases added to the lysate by the PreON. Digestion reaction is stopped.
3. Each lysate protease mix is transferred to a spin column-based spin column-based cartridge in a rotor adapter.
4. Peptides bind to the membrane. During washing steps of the PreOmics spin column-based spin column-based cartridge salts and contaminants are removed.
5. The spin column-based cartridge is transferred to a microcentrifuge tube for elution of purified peptides.


For an example workflow using a PreOmics iST Kit, see the flowchart on the next page.

3.1 External features of the PreON

Front view of the PreON

Pulled-out touchscreen

Rear view of PreON

Rear view of PreON

- | | | | |
|---|--------------------|---|--|
| 1 | Touchscreen | 5 | 2x USB ports on the left side of the touchscreen;
2x USB ports behind the touchscreen
(Wi-Fi module plugged into 1 USB port) |
| 2 | Hood | 6 | RJ-45 Ethernet port |
| 3 | Waste drawer | 7 | Power cord socket |
| 4 | Front power switch | 8 | Cooling air outlet |

Touchscreen

The PreON is controlled using a swivel-mounted touchscreen. The touchscreen allows the user to select and run protocols, install or delete protocols, save data files and operate the shaker and centrifuge individually. Information displayed on the touchscreen guides the user through worktable setup after a protocol has been selected. During sample processing, the touchscreen shows the protocol status and remaining time.

Pulled-out touchscreen

Hood

The PreON hood protects users from the moving robotic arm and from potentially infectious material placed on the worktable. The hood can be manually opened to gain access to the worktable. During operation of the PreON, the hood must remain closed and should only be opened when instructed to do so by the software. Therefore, a hood lock has been implemented to avoid inappropriate opening.

WARNING

Moving parts

[W18]

To avoid contact with moving parts during operation of the PreON, the instrument must be operated with the hood closed.

If the hood sensor or lock is not functioning correctly, contact PreOmics Technical Services.

Power switch

The power switch is located at the front right of the PreON. To power ON the PreON, press the switch. The beeper sounds, and the startup screen appears. The instrument then automatically performs initialization tests.

To conserve energy, the PreON can be powered OFF when not in use. To power OFF the PreON, press the switch.

Note: After powering OFF the PreON, you must wait for a few seconds before you power ON the instrument again. The system might fail to start if you do not allow the PreON to rest for a few seconds before powering ON.

RJ-45 Ethernet port

The RJ-45 Ethernet port located at the back of the instrument beside the power cord socket, is only used to connect the PreON to local area network via cable.

USB ports

The PreON has four USB ports. Two are located on the left of the touchscreen and two are located behind the touchscreen.

The USB ports located on the left of the touchscreen allow connection of the PreON to a USB stick. Data files, such as support package, protocols or report files, can be transferred via the USB port from the PreON to the USB stick. The USB ports can also be used to plug in an external 2D bar code reader.

New protocols and software versions, available at www.preomics.com/downloads, can be downloaded to the USB stick and transferred to the PreON via the USB port.

Important: The USB port is only for use B devices provided by PreOmics. Do not connect other devices to this port.

Important: Do not remove the USB stick while downloading or transferring data or software to or from the instrument.

Waste drawer

Used disposable filter-tips are discarded through two slots in the worktable and collected in the waste drawer.

CAUTION

Damage to the instrument

[C3]

Empty the tip disposal container prior to use to prevent a tip jam in the waste drawer. Failure to empty the waste container may block the robotic arm that could cause run failure or instrument damage.

WARNING

Hazardous chemicals and infectious agents

[W16]

The waste may contain toxic material and must be disposed of properly. Refer to your local safety regulations for proper disposal procedures.

WARNING

Hazardous chemicals

[W13]

Some chemicals used with this instrument may be hazardous or may become hazardous after completion of the protocol run.

Always wear safety glasses, gloves, and a lab coat.

The responsible body (e.g., laboratory manager) must take the necessary precautions to ensure that the surrounding workplace is safe and that the instrument operators are not exposed to hazardous levels of toxic substances (chemical or biological) as defined in the applicable Safety Data Sheets (SDSs) or OSHA³, ACGIH⁴ documents.

Venting for fumes and disposal of wastes must be in accordance with all national, state, and local health and safety regulations and laws.

³ OSHA : Occupational Safety and Health Administration (United States of America).

⁴ ACGIH : American Conference of Government Industrial Hygienists (United States of America).

WARNING

Samples containing infectious agents

[W12]

Some samples used with this instrument may contain infectious agents. Handle such samples with the greatest of care and in accordance with the required safety regulations.

Always wear safety glasses, 2 pairs of gloves, and a lab coat.

The responsible body (e.g., laboratory manager) must take the necessary precautions to ensure that the surrounding workplace is safe, and that the instrument operators are suitably trained and not exposed to hazardous levels of infectious agents as defined in the applicable Safety Data Sheets (SDSs) or OSHA⁶, ACGIH⁵ documents.

Venting for fumes and disposal of wastes must be in accordance with all national, state, and local health and safety regulations and laws.

Power cord socket

The power cord socket is located at the rear right of the PreON and allows connection of the PreON to a power outlet via the supplied power cord.

WARNING

Electrical hazard

[W8]

Any interruption of the protective conductor (earth/ground lead) inside or outside the instrument or disconnection of the protective conductor terminal is likely to make the instrument dangerous.

Intentional interruption is prohibited.

Lethal voltages inside the instrument

When the instrument is connected to line power, terminals may be live and opening covers or removing parts is likely to expose live parts.

WARNING

Damage to electronics

[W9]

Before powering ON the instrument, make sure that the correct supply voltage is used.

Use of incorrect supply voltage may damage the electronics.

To check the recommended supply voltage, refer to the specifications indicated in the type plate of the instrument.

WARNING

Risk of electric shock

[W10]

Do not open any panels on the PreON.

Risk of personal injury and material damage

Only perform maintenance that is specifically described in this user manual.

⁶ OSHA : Occupational Safety and Health Administration (United States of America).

⁵ ACGIH : American Conference of Government Industrial Hygienists (United States of America).

Cooling air outlet

Cooling air outlets are located at the rear left side of the PreON and allow cooling of the internal components of the PreON.

CAUTION


Risk of overheating

[C11]

To ensure proper ventilation, maintain a minimum clearance of 10 cm at the sides and rear of the PreON.

Slits and openings that ensure the ventilation of the PreON must not be covered.

3.2 Internal features of the PreON

Internal view of the PreON.

- | | | | |
|---|--------------------------|---|---|
| 1 | Centrifuge lid | 6 | Microcentrifuge tube slots |
| 2 | Centrifuge | 7 | Slots for tip racks |
| 3 | Shaker | 8 | Disposal slots for tips and columns |
| 4 | Liquid waste | 9 | Robotic arm (includes 1 channel pipettor, gripper, ultrasonic and optical sensor) |
| 5 | Tip sensor and hood lock | | |

Centrifuge

The centrifuge is equipped with 12 swing-out buckets, each of which can hold a disposable rotor adapter. Up to 12 samples can be processed per run. For ease of use and high process safety, a grey line marks the side of the bucket that must face toward the center of the rotor. All centrifuge buckets must be mounted before starting a run, regardless of the number of samples to be processed. Be sure to follow the loading instructions provided by the software to ensure correct loading of the centrifuge.

The centrifuge can also be operated individually via the touchscreen.

Note: Be sure to follow the centrifuge loading instructions provided by the software.

WARNING

Risk of personal injury and material damage

[W5]

To prevent plastic crash, load the tubes properly. After a plastic crash, sharp plastic particles, toxic and aggressive reagents could be inside the centrifuge. Be careful when handling items inside the centrifuge and thoroughly clean the centrifuge according to the maintenance procedure.

WARNING


Risk of personal injury and material damage

[W20]

Raise the centrifuge lid carefully. The lid is heavy and may cause injury if it falls.

Rotor adapter

A disposable rotor adapter holds a PreOmics spin column-based cartridge placed in a spin column-based cartridge adaptor and a microcentrifuge tube in a centrifuge bucket during sample processing. For ease of use and high process safety, the rotor adapters are designed so that they fit into a centrifuge bucket only in the correct orientation. Spin column-based cartridge and microcentrifuge tube lids are held securely in slots at the edge of the rotor adapter.

Microcentrifuge tube position

Middle position

Wash position

Slot for spin column-based cartridge lid

Assembly of a rotor adapter

The wash position of the rotor adapter is open at the bottom, enabling wash buffers to flow through and collect at the bottom of the rotor adapter during centrifugation. The other two positions in the rotor adapter are closed. Be sure to follow the loading instructions provided by the software.

WARNING

Risk of personal injury and material damage

[W4]

Do not use damaged rotor adapters. The rotor adapters can only be used once. High g forces exerted in the centrifuge can cause damage to used rotor adapters.

WARNING

Risk of personal injury and material damage

[W5]

To prevent plastic crash, load the tubes properly. After a plastic crash, sharp plastic particles, toxic and aggressive reagents could be inside the centrifuge. Be careful when handling items inside the centrifuge and thoroughly clean the centrifuge according to the maintenance procedure.

Shaker

The heated orbital shaker enables fully automated lysis of up to 12 or 16 samples. Sample tubes are placed into a rack that fits onto the shaker adapter. The lid of each microcentrifuge tube is held in a slot at the edge of the shaker rack. This ensures that microcentrifuge tubes cannot be displaced during sample processing and allows shaker loading to be checked. The positions on the shaker rack are numbered to enable easy loading.

Shaker rack for 12 samples

Sample tube lids are held in slots at the edge of the shaker rack

Shaker rack for 16 samples

Sample tube lids are held in slots at the edge of the shaker rack

Note: For shaker loading follow instructions provided by the software.

WARNING

Risk of personal injury and material damage

[W5]

To prevent plastic crash, load the tubes properly. After a plastic crash, sharp plastic particles could be inside the centrifuge. Be careful when handling items inside the centrifuge.

WARNING

Hot surface

[W21]


The shaker can reach temperatures of up to 83°C (181.4°F). Avoid touching it when it is hot. Carefully remove the samples after a run.

Additional reagent racks - Reagent work decks for position Insert 3

Insert 3 for screwcap microtubes:

The additional reagent work deck for position Insert 3 is able to host 15+1 screwcap microtubes. Liquid is aspirated from the microtubes by the pipetting system.

Note: Screw-cap microtubes are used with the PreON.

Insert 3 with 16 microtubes positions and waste bottle

Insert 3 for 16 plex TMT:

The additional reagent work deck for position Insert 3 is able to host 16 TMT tubes and 4 screwcap microtubes. Liquid is aspirated from the microtubes by the pipetting system.

Note: Screw-cap microtubes are used with the PreON.

Insert 3 with 16 microtubes positions

WARNING

Risk of fire or explosion

[W6]

When using flammable liquids on the PreON (Acetonitrile, Isopropyl alcohol), handle such liquids carefully and in accordance with the required safety regulations. If liquid has been spilled, wipe it off and leave the PreON hood open to allow flammable vapors to disperse.

WARNING

Risk of explosion

[W7]

The PreON is intended for use with reagents and substances supplied with PreOmics kits or other than outlined in respective Information for use. Use of other reagents and substances may lead to fire or explosion.

WARNING

Hot surface

[W21]

The shaker can reach temperatures of up to 83°C (181.4°F). Avoid touching it when it is hot. Carefully remove the samples after a run.


Tip sensor

During sample preparation, the tip sensor checks that the tip adapter has picked up a tip.

Tip rack slots

Two tip racks can be placed on the PreON worktable. Tips can be purchased in prefilled tip racks holding 200 µl filter-tips.

Note: Only filter-tips designed for use with the PreON must be used.

Tip disposal slots

Disposable filter-tips are alternately discarded through each of the tip disposal slots into the waste drawer. This prevents discarded tips from piling up in the waste drawer.

Robotic arm

The robotic arm provides accurate and precise positioning of the robotic gripper and pipetting system on the PreON worktable and includes an optical and ultrasonic sensor.

WARNING

Moving parts

[W18]

To avoid contact with moving parts during operation of the PreON, the instrument must be operated with the hood closed.

If the hood sensor or lock is not functioning correctly, contact PreOmics Technical Services.

WARNING


Moving parts

[W17]

Avoid contact with moving parts during operation of the PreON. Under no circumstances should you place your hands under the robotic arm when it is lowering. Do not attempt to move any tip racks or tubes whilst the instrument is operating.

Robotic gripper

The robotic gripper transfers spin column-based cartridges. During transfer of a spin column-based cartridge, a stabilizing rod holds the rotor adapter in place, ensuring it remains properly seated in the centrifuge bucket. The robotic gripper is behind the panel covering the robotic arm.

Robotic gripper

Stabilizing rod

Fully automated spin column-based cartridge processing

Pipetting system

The PreON is equipped with a single-channel pipetting system that moves in the X, Y and Z directions. The dilutor, fitted with a tip adapter, is connected to a precision syringe pump, which enables accurate transfer of liquids. The tip adapter allows aspiration and dispensing of liquid through an attached disposable tip. Disposable filter-tips (200 µl) are used for sample processing to minimize the risk of cross-contamination.

Optical sensor

During the load check, the optical sensor checks that the number of rotor adapters corresponds to the number of samples in the shaker and that both shaker and rotor are correctly loaded.

Ultrasonic sensor

During the load check, the ultrasonic sensor checks that the reagent microtubes contain enough liquid for the protocol run. The ultrasonic sensor also checks the type of tips loaded on the worktable and whether there are enough tips for the protocol run.

Note: The ultrasonic sensor has a black beam columnator. If, for any reason, this beam columnator falls off or is missing, the instrument will display an error message to inform the user that the beam columnator is missing and that runs cannot be started. To replace the beam columnator, it needs to be manually adjusted to its original position (see image below).

Black beam columnator (see red circle) of the ultrasonic sensor

External hand scanner

The PreON is equipped with a 2D handheld scanner model Gryphon GD4400 to allow kit bar code and sample bar code scanning. Appendix A – Technical data on page 112

WARNING

Risk of personal injury

Hazard Level 2 laser light: Do not stare into the light beam when using handheld bar code scanner.

[W28]

4 Installation Procedures

This section provides instructions on unpacking, packing and installing the PreON.

4.1 Installation environment

4.1.1 Site requirements

The PreON must be located out of direct sunlight, away from heat sources and away from sources of vibration and electrical interference. Refer to Appendix A – Technical data on page 112 for the operating conditions (temperature and humidity). The site of installation should be free of excessive drafts, excessive moisture and excessive dust and should not be subject to large temperature fluctuations.

The PreON must be used only in well ventilated areas.

Use a level workbench that is large enough and strong enough to accommodate the PreON. Refer to Appendix A – Technical data for the weight and dimensions of the PreON.

Ensure that the workbench is dry, clean and vibration-proof and has additional space for accessories.

The PreON must be placed within approximately 1.5 m of a properly grounded (earthed) AC power outlet. The power line to the instrument should be voltage regulated and surge protected. Ensure that the PreON is positioned so that it is easy to access the power connector at the back of the instrument and the power switch on the front at all times, and that it is easy to power the instrument OFF and disconnect it.

Note: It is recommended to plug the instrument directly into its own power outlet and not to share the power outlet with another lab equipment. Do not place the PreON on a vibrating surface or near vibrating objects.

WARNING

Explosive atmosphere

[W11]

The PreON is not designed for use in an explosive atmosphere.

WARNING

Vapors from volatile solvents

[W14]

Some of the solvents used with the PreON are volatile. Use the PreON only with adequate ventilation according to the applicable Safety Data Sheets (SDSs) or OSHA,* ACGIH† documents.

CAUTION

Risk of overheating

[C11]

To ensure proper ventilation, maintain a minimum clearance of 10 cm at the sides and rear of the PreON.

Slits and openings that ensure the ventilation of the PreON must not be covered.

WARNING

Risk of personal injury and material damage

[W2]

The PreON is too heavy to be lifted by one person. To avoid personal injury or damage to the instrument, do not lift the instrument alone.

CAUTION

Damage to the instrument

[C6]

Direct sunlight may bleach parts of the instrument and cause damage to plastic parts.

The PreON must be located out of direct sunlight.

4.1.2 Power requirements

The PreON operates at: 100–240 V AC, 50/60 Hz, 650 VA

Ensure that the voltage rating of the PreON is compatible with the AC voltage available at the installation site. Mains supply voltage fluctuations are not to exceed 10% of nominal supply voltages.

WARNING

Damage to electronics

[W9]

Before powering ON the instrument, make sure that the correct supply voltage is used.

Use of incorrect supply voltage may damage the electronics.

To check the recommended supply voltage, refer to the specifications indicated in the type plate of the instrument.

WARNING

Electrical hazard

[W8]

Any interruption of the protective conductor (earth/ground lead) inside or outside the instrument or disconnection of the protective conductor terminal is likely to make the instrument dangerous.

Intentional interruption is prohibited.

Lethal voltages inside the instrument

When the instrument is connected to line power, terminals may be live and opening covers or removing parts is likely to expose live parts.

4.1.3 Grounding requirements

To protect operating personnel, the PreON shall be correctly grounded (earthed). The instrument is equipped with a 3-conductor AC power cord that, when connected to an appropriate AC power outlet, grounds (earths) the instrument. To preserve this protection feature, do not operate the instrument from an AC power outlet that has no ground (earth) connection.

WARNING

Electrical hazard

[W8]

Any interruption of the protective conductor (earth/ground lead) inside or outside the instrument or disconnection of the protective conductor terminal is likely to make the instrument dangerous.

Intentional interruption is prohibited.

Lethal voltages inside the instrument

When the instrument is connected to line power, terminals may be live and opening covers or removing parts is likely to expose live parts.

4.2 Installation of AC power cord

The AC power cord connects to the socket located at the rear of the PreON, and the other end to the AC power outlet.

4.3 Unpacking the PreON

WARNING

Risk of personal injury and material damage

[W2]

The PreON is too heavy to be lifted by one person. To avoid personal injury or damage to the instrument, do not lift the instrument alone.

Important: Only qualified service personnel is allowed to unpack, install the PreON & re-pack the device for returns.

14. Before unpacking the PreON, move the package to the site of installation and check that the arrows on the package point upward. In addition, check whether the package is damaged. In case of damage, contact PreOmics Technical Service.

15. Open the top of the transportation box to remove the PreON power cord before lifting the box.

16. Remove the black foam protector lid and lift the box.

17. When lifting the PreON, slide your fingers under both sides of the workstation and keep your back straight.

Important: Do not hold the touchscreen display while unpacking or lifting the PreON, this might damage the instrument.

18. Check if the packing list document is included after unpacking the PreON.

19. Read the packing list to check that you have received all items. If anything is missing, contact PreOmics Technical Services.

20. Check that the PreON is not damaged and that there are no loose parts. If anything is damaged, contact PreOmics Technical Services. Make sure that the PreON has equilibrated to ambient temperature before operating it.

21. Retain the package in case you need to transport the PreON in the future. Refer to section 4.6, Packing the PreON for more details. Using the original package minimizes the possibility of damage during transportation of the PreON.

4.4 Installing the PreON

Important: Only qualified service personnel is allowed to unpack, install the PreON & re-pack the device for returns.

This section describes important actions that must be performed before operating the PreON. These actions include:

- Removal of the PreON accessories.
- Removal of the protective film from the PreON hood.
- Removal of the protector for the robotic arm.
- Connection of the power cord to the back of the PreON.
- Powering ON the PreON.
- Installation of the centrifuge rotor and buckets.

4.4.1 Remove PreON accessories

Remove the power cord and all accessories from the foam packing material on top of the PreON.

Remove the USB stick, rotor key, rotor nut, Allen key and shaker rack plugs from the waste drawer.

4.4.2 Remove the foam protector

During transportation, a foam protector prevents the movable parts of the PreON from moving along the X- and Y-axes. Before using the instrument, this protector must be removed.

Foam protector for robotic arm.

To remove the foam protector, gently pull the foam protector towards you (see picture above). After removing the protector for the robotic arm, make sure to close the PreON hood.

4.4.3 Installation of AC power cord

1. Remove the power cord from the foam packing material on top of the PreON.
Note: Only use the power cord that is supplied with the PreON.
2. Ensure that the power switch is set to OFF: outer position is OFF and inner position is ON.
3. Check that the voltage rating on the label at the back of the PreON matches the voltage available at the installation site.

4. Plug the power cord into the instrument power-cord socket.
5. Plug the power cord into a grounded power outlet.

WARNING

Damage to electronics

[W9]

Before powering ON the instrument, make sure that the correct supply voltage is used.

Use of incorrect supply voltage may damage the electronics.

To check the recommended supply voltage, refer to the specifications indicated in the type plate of the instrument.

WARNING

Electrical hazard

[W8]

Any interruption of the protective conductor (earth/ground lead) inside or outside the instrument or disconnection of the protective conductor terminal is likely to make the instrument dangerous.

Intentional interruption is prohibited.

Lethal voltages inside the instrument

When the instrument is connected to line power, terminals may be live and opening covers or removing parts is likely to expose live parts.

4.4.4 Powering ON the PreON

1. Check that the PreON operates properly:
2. Make sure that the hood of the PreON is closed.
3. Power ON the PreON using the power switch (inner position is ON; outer position is OFF).

The beeper sounds and the startup screen appears.

The instrument automatically performs initialization tests.

The centrifuge lid opens.

4. If there is an initialization error, check that the power cord is properly connected to the PreON and power outlet. Retry the initialization process. If the problem persists, contact PreOmics Technical Services.

Note: If you powered ON the PreON before closing the hood, the instrument will not perform an initialization and the centrifuge lid will not open. Power OFF the PreON, close the hood, and then power ON the PreON again.

4.4.5 Installing the centrifuge rotor and buckets

The centrifuge rotor and buckets are pre-installed in the PreON. When setting up the PreON for the first time, turn on the instrument and remove the transportation foam inserts from the centrifuge after it opens. In case you have manually removed the centrifuge rotor and buckets (e.g., during maintenance), follow the instructions below to re-install them.

The rotor can be mounted in only one orientation: a pin on the rotor shaft fits into a notch on the underside of the rotor directly underneath rotor position 1. Line up position 1 of the rotor with the pin on the rotor shaft and carefully lower the rotor onto the shaft. Install the rotor nut on top of the rotor and tighten using the rotor key supplied with the PreON. Make sure that the rotor is securely seated. Check that all buckets are properly suspended and can swing freely.

Note: The side of the rotor bucket that must face towards the rotor shaft is marked with a grey line to help prevent buckets from being loaded incorrectly.

Important: All centrifuge buckets must be mounted before starting a protocol run, even if fewer than 12 samples are to be processed.

WARNING

Risk of personal injury and material damage

[W25]

To prevent the rotor nuts from loosening during operation of the centrifuge, securely tightened the nuts using the rotor key supplied with PreON.

WARNING

Risk of personal injury and material damage

[W20]

Raise the centrifuge lid carefully. The lid is heavy and may cause injury if it falls.

CAUTION

Damage to the instrument

[C7]

The PreON must not be used if the centrifuge lid is broken, or if the lid lock is damaged.

Make sure that no loose material is inside the centrifuge during operation.

Make sure that the rotor is installed correctly and that all buckets are properly mounted, regardless of the number of samples to be processed. Load the rotor only as instructed by the software.

Only use rotors, buckets, and consumables designed for use with the PreON. Damage caused by use of other consumables will void your warranty.

We recommend replacing the centrifuge after 20,000 cycles, which is equivalent to 9 years of usage with two runs per day for 220 days each year. For more information contact PreOmics Technical Support.

4.5 Configuration of the PreON

When using the PreON for the first time, it is recommended to define the required settings. Other settings can be made later when needed.


For details on using the touchscreen and software, refer to section 5.1 PreON software on page 53

To configure the PreON, follow the steps below.

1. Close the instrument hood.
2. Press the power switch to the inner position to power ON the instrument. The startup screen appears and the beeper sounds (if enabled in the sound settings). The instrument automatically performs the initialization tests. If the centrifuge lid is closed, it will open.
3. Initially, only one user account is available: the pre-installed default user. Press **OK** on the touchscreen to confirm the message.
4. Initially, if no maintenance has been recorded yet, the maintenance status is initialized using default values. Press **OK** to confirm the message.

5. The **Login** screen appears after the initialization.

PREOMICS

Please log in to start


User ID

Password

Log in

Login screen

6. Initially, only one default user is available. In this case, enter '**Admin**' in the **User ID** and 'admin' in the **Password** field using the on-screen keyboard. Touch in the entry field to open the on-screen keyboard. For details, refer to section 5.1.1 Entering text and numbers.
7. Press **Log in**. The **Home** (main) screen appears.

Home (main) screen

8. If you need to return to the **Home** screen from another screen, press the **Setup** icon.

4.5.1 Setting basic system data

This section describes how to set the following:

- Name for the PreON
- Current date and time

When using the PreON for the first time, it is recommended to set the current date and time.

1. On the main screen, press the **Configuration** icon.

2. Press the **System** tab. This tab is only available for users assigned the Administrator role.

PREOMICS

Setup Tools **Configuration** Logout

System Process Users Protocols Settings LAN Shaker

Device name
PreON

Date May 02, 2022 Time 08 : 37

Serial number 0
Software version ?
For the latest software version, please check PREOMICS.COM. Use the USB stick provided to transfer the downloaded update on the instrument.

USB Connected

Reset Rotor Counter Update Software

May 02, 2022, 08:37 User: Admin Admin

3. Optionally, enter a name for the PreON. The device name serves as the network/host name when connecting the instrument to the network.

The name can have up to 24 characters: letters A–Z, a–z, digits 0–9 and a hyphen (-).

The name must start with a letter and cannot end with a hyphen (-).

4. In the **Date** and **Time** fields, select the current date and enter the current time for the instrument. These are used to track the start and end time of a run and are also part of the run report. Date and time are not synchronized using the network. To change the date, touch the calendar icon (📅) and select the date.
5. Use the left and right arrow icons to change the month then touch the current day and press **OK** to confirm.

PREOMICS

System Process Users Protocols Settings LAN Shaker

Device name
PreON

Date May 02, 2022

Serial number 0
Software version ?
For the latest software version, please check PREOMICS.COM. Use the USB stick provided to transfer the downloaded update on the instrument.

USB Connected

Reset Rotor Counter Update Software

May 02, 2022, 08:37 User: Admin Admin

May 2022

Mon	Tue	Wed	Thu	Fri	Sat	Su
25	26	27	28	29	30	1
2	3	4	5	6	7	8
9	10	11	12	13	14	15
16	17	18	19	20	21	22
23	24	25	26	27	28	29
30	31	1	2	3	4	5

Cancel OK

4.5.2 User settings

This section describes the optional settings that can be defined by each user:

- Audio volume
- Display brightness
- Internal light

The settings apply to the current user.

1. On the main screen, press the **Configuration** icon.

2. Press the **Settings** tab. This tab is available for all users.

3. To adjust the audio volume or display brightness, touch the desired position on the virtual slider on the screen. For audio volume, a sound is played with the set volume.
4. Check the box next to **Use internal light** to turn on the LED light inside the instrument. Uncheck the box to turn it off.
5. To return to the main screen, press the **Setup** icon.

4.5.3 Managing users

The PreON is provided with a User Management feature. This feature enables you to set up multiple users with two different roles: administrator and operator. When you use the PreON for the first time, a default user named Admin is already pre-installed and configured with both roles assigned. The user management feature is only available for users assigned the Administrator role.

Setting up a new user

1. Press the **Configuration** icon on the main screen.

2. Press the **Users** tab.

The configured users are shown in the table. Each row contains the data for one user.

The screenshot shows the PREOMICS Configuration menu with the 'Users' tab selected. Below the menu is a table of configured users.

User Id	First Name	Last Name	Role(s)	Edit / Delete
Admin	Admin	Admin	Administrator, Operator	

Below the table, there is a checkbox labeled 'Show only activated user profile' which is checked, and a 'New ...' button.

May 02, 2022, 08:40 User: Admin Admin

Note: It is recommended to create at least one other user with the administrator role.

3. Press **New** to add a new user.
4. Enter the respective data for the new user. Keep the **Activate User** box checked.

The screenshot shows the 'Add User' dialog box with the following fields and options:

- User Id:
- First name:
- Last name:
- E-mail:
- Enter password:
- Confirm password:
- Role selection: ☐ Administrator, ☐ Operator
- Activate User: ☒
- Buttons: Cancel, OK

The **User ID**, **First name** and **Last name** fields are mandatory. These fields may contain up to 30 letters and numerical characters. The user ID must be unique for each user profile. It must contain at least one letter and cannot contain blank spaces. The user ID is used for logging in and is printed on run reports. The first and last name are displayed on the touchscreen for the currently logged in user.

The **Password** field is mandatory and must contain 8–40 letters or numerical characters. Enter the same password into the **Confirm password** field.

Select the user role: **Administrator** and/or **Operator**. The operator may only use the instrument, while the administrator is also allowed to configure the system. One user can have both roles assigned at the same time. The default user **Admin** has both user roles assigned.

The **E-mail** address field is optional. The system does not confirm whether the e-mail address entered is valid.

5. Press **OK** to save the new user.

Changing data for an existing user

1. Press the **Configuration** icon on the main screen.

2. Press the **Users** tab.

The configured users are shown in the table. Each row contains the data for one user.

User Id	First Name	Last Name	Role(s)	Edit / Delete
Admin	Admin	Admin	Administrator, Operator	
jane	Jane	Doe	Administrator	

☒ Show only activated user profile New ...

May 02, 2022, 08:42 User: Admin Admin

3. In the user profile row, press the **Edit** () icon.
4. A screen will appear showing the current information of the user. Edit the information as necessary.

Edit User

Anonymous user ID: First name: Last name:

E-mail:

Enter password: ☐ Administrator ☒ Operator

Confirm password:

☒ Activate User Cancel OK

The user's password will not be displayed. If you touch the password field, the existing password will be cleared, and a new password must be entered and confirmed.

5. To confirm the changes, press **OK**. To close the dialog and discard the changes, press **Cancel**.

Deleting or temporarily deactivating a user

1. To delete a user, press the **Delete** icon () in the user profile row. It is not possible to delete the currently logged-in administrator.
2. To temporarily deactivate a user, press the **Edit** icon () in the user profile row. Uncheck the **Activate User** box. It is not possible to deactivate the currently logged-in administrator.

3. To reactivate a user profile, press the **Edit** icon () in the user profile row. Check the **Activate user** box.

Note: If a user tries to login with the wrong password, the user profile will be automatically deactivated after 10 failed login attempts.

Deactivated users can be displayed in the users list by unchecking the box **Show only activated user profiles**. If the box is unchecked, all user profiles will be listed.

If the pre-installed administrator user **Admin** has more than 10 failed login attempts, the password can be reset. To do so, press **Reset admin user**. Then log in again with the default user ID **Admin** and password **admin**.

PREOMICS

Please log in to start

User ID

Password

 You have entered wrong credentials.

Log in

4.5.4 Changing password

Every user can change their own password.


A user with the role **Administrator** is allowed to change the password for every user by editing the user profile. Refer to section Changing data for an existing user on page 42 for more details. Passwords are never displayed in this process, so the administrator cannot view a user's password.

Changing your password

1. Press the **Configuration** icon on the main screen.

- For users with the role **Operator**, the **Password** tab is automatically active.

- Enter the old password into the **Old password** field. Touch the field to open the on-screen keyboard.
Enter a new password into the **New password** field and re-enter the new password in the **Confirm new password** field. The **Password** must contain 8–40 letters or numerical characters.
- Press **OK** to save the new password. Press **Cancel** to discard any changes and to keep the old password.
- To return to the main screen, press the **Setup** icon.

4.5.5 Managing protocols

Commonly used PreOmics standard protocols are installed on the PreON upon delivery. The range of PreOmics standard protocols is continuously expanding. To learn more, please reach out to your local PreOmics contact. If unsure, get in touch at preomics.com/contact. PreOmics's Application Lab Specialists can also customize these protocols or develop new protocols depending on your needs. Protocols that are no longer required can be removed from the PreON.


Protocols can only be managed by users assigned the Administrator role.

Installing new protocols

- On a computer running Microsoft® Windows®, download the new protocols provided to you by PreOmics. Use the USB stick that was shipped with the PreON to transfer the protocol files to the instrument. Copy the new protocol zip file(s) to the USB stick into the **Protocol_Upload** directory. Do not unzip the files. Make sure to use the correct directory, otherwise, the PreON will not find the protocols.
Note: Do not rename or modify the protocol files. Otherwise, they cannot be used.
- Connect the USB stick to the PreON using one of the USB ports at the left of the touchscreen.
- Select the **Configuration** icon.

4. Press the **Protocols** tab.

5. Press **Copy from USB**.

6. A message is displayed indicating how many protocols are found on the USB stick. Press **Yes** to start uploading.

All protocol zip file(s) in the **New Protocols** folder will be installed.

Note: Already installed protocols will not be overwritten. If you attempt to re-install an existing protocol, a message will appear indicating that not all protocols could be copied.

7. Wait until the transfer is complete. A message is displayed when the transfer is complete.
8. Remove the USB stick and power OFF the PreON.
9. Wait a few seconds and then power ON the PreON. To use the new protocols, log in again.
- To transfer all installed protocols to the USB stick, press **Copy to USB**.

Deleting protocols

Important: Before deleting, you should back up the protocols on the USB stick provided with the instrument.

1. Select the **Configuration** icon.

2. Press the **Protocols** tab.

3. To delete all protocols installed on the instrument, press **Delete All**. It is not possible to delete a single protocol from the PreON.

Note: After deleting all protocols, a selective upload of protocols will reduce the options during run setup.

4.5.6 Updating software

Please reach out to your local PreOmics contact or get in touch at preomics.com/contact to learn if an updated software version is available for download.

The software can only be updated by users assigned the role Administrator. It is recommended to download all run reports before updating the software, because run reports and support packages will be lost during software update (see section 5.7 Saving run reports to the USB stick).

1. On the main screen, press the **Configuration** icon.

2. Press the **System** tab.
3. The currently installed software version is shown at the right.

4. On a computer running Microsoft Windows, download the software to the main folder of the USB stick provided with the PreON. Wait until the download has completed before removing the USB stick.

Note: After download, make sure the following 3 files are in the main folder of the USB stick:

preon1.bin
preon2.bin
preon-<version>.tar.gz

The update will not work if one of the files is missing or has been renamed. Make sure that only the files for one software version are in the root folder of the USB stick.

5. Connect the USB stick to the instrument using one of the USB ports at the left of the touchscreen.

Important: Make sure that all run reports and support packages have been backed up before proceeding to the next step.

- Press **Update Software** to start the software update. Follow the instructions on the screen.

PREOMICS

Setup Tools Configuration Logout

System Process Users Protocols Settings LAN Shaker

Device name
PreON

Date
May 02, 2022

Time
13 : 07

Serial number
0

Software version
1.10.7

For the latest software version, please check PREOMICS.COM. Use the USB stick provided to transfer the downloaded update on the instrument.

USB
Connected

Reset Rotor Counter

Update Software

May 02, 2022, 13:07 User: Admin Admin

- The software update tool is displayed. Press **Start update** to start the software update.

PREOMICS Preon Software Update Tool 1.0.0

Update Preon Software

i All internal data and settings including run reports and support packages get lost. To backup these files unplug the USB stick, restart the device and download the log files and run reports in the "Data" tab. Press start to update the Preon software

Replacing Version 1.0.0 with 1.1.0

0%


Start Update Cancel

▼

Press **Cancel** if you do not want to update the software. In this case, the instrument will initialize without updating the software.

- Wait until the update has completed.

- Press the arrow button to open the **Details** screen.

- Compare the MD5 checksum on the screen with the checksum provided on the software download page. If the checksums are not identical, contact PreOmics Technical Services.
- Press **Reboot** to continue. The instrument will initialize with the updated software.
- When instructed by the screen, remove the USB stick from the USB port.
- Use a computer running Microsoft Windows to delete the previously downloaded software files from the USB stick.

4.6 Packing the PreON

Warning: The PreON is too heavy to be lifted by one person. To avoid personal injury or damage to the instrument, do not lift the instrument alone.

WARNING

Risk of personal injury and material damage

[W2]

The PreON is too heavy to be lifted by one person. To avoid personal injury or damage to the instrument, do not lift the instrument alone.

Important: Only qualified service personnel is allowed to unpack, install the PreON & re-pack the device for returns.

Before transporting the PreON, the instrument must first be decontaminated. Refer to section 6.7, Decontaminating the PreON for more details. Then prepare the instrument as follows:

- Prepare the packing material. Materials required are the cardboard carton, the pallet with foam blocks, the foam lid and the foam protector for the robotic arm.
Note: The centrifuge lid must be open to allow access to the inside of the centrifuge. If the lid is not open, perform steps 2–5 below. If the lid is already open, proceed to step 6.
- Close the instrument hood.
- On the **Main** screen, press the **Tools** button.
- In the Tools menu, press the Run Modules tab.
- In the Run Modules menu, press the Centrifuge tab and then press Open Centrifuge Lid.
- Undo the rotor nut on top of the rotor using the rotor key, and carefully lift the rotor off the rotor shaft.

7. Close the hood.
8. On the **Main** screen, press the **Tools** button.
9. In the **Tools** menu, press the **Run Modules** tab.
10. In the **Run Modules** menu, press the **Centrifuge** tab and then press **Close Centrifuge Lid**.
11. When the centrifuge lid is closed, power OFF the PreON and open the hood.
12. Insert the foam protector into the front of the instrument.
13. Press the foam down between the centrifuge and the robotic arm.

Foam protector inserted between the centrifuge and the robotic arm

14. Push the foam until the rear end touches the back wall of the instrument. Ensure that the arm is held firmly in place and cannot move.
15. Make sure that the PreON hood can be closed properly. The hood should lightly brush against the foam.
16. Place the accessories into the waste drawer. The following accessories should be packed in air cushion bags:
 - Rotor key
 - Allen key
 - Rotor nut
 - USB stick
 - Shaker rack plugs
17. Place the handheld scanner into dedicated handheld scanner box.
18. Place the PreON onto the pallet and put the black foam lid over the top of the instrument. Place the box onto the instrument.

Important: When lifting the PreON, slide your fingers under both sides of the instrument and keep your back straight.

Important: Do not hold the touchscreen display while lifting the PreON, as this might damage the instrument.

WARNING

Risk of personal injury and material damage

[W2]

The PreON is too heavy to be lifted by one person. To avoid personal injury or damage to the instrument, do not lift the instrument alone.

19. Place the accessories into the black foam lid. The following accessories should be packed in air cushion bags:
 - Rotor with swing-out buckets
 - Power cord
20. Seal the outside edges of the carton with tape to protect against moisture.

Note: Using the original package minimizes potential damage during transportation of the PreON.

5 Operating Procedures

This section describes how to operate the PreON.

Before proceeding, we recommend that you familiarize yourself with the features of the instrument by referring to sections 3.1 and 3.2.

Important: The PreON is designated for use with PreOmics spin column-based cartridges only. Use of other types of spin column-based cartridges or chemistries will void your warranty.

CAUTION

Damage to the instrument

[C2]

Only use PreOmics spin column-based cartridges and PreON specific consumables with the PreON. Damage caused by use of other types of spin column-based cartridges or chemistries will void your warranty.

The hood of the PreON must remain closed during operation of the instrument. Only open the hood when instructed to do so by the software.

WARNING

Moving parts

[W18]

To avoid contact with moving parts during operation of the PreON, the instrument must be operated with the hood closed.

If the hood sensor or lock is not functioning correctly, contact PreOmics Technical Services.

WARNING

Risk of personal injury and material damage

[W3]

Do not attempt to move the PreON during operation.

WARNING

Risk of personal injury and material damage

[W1]

Improper use of the PreON may cause personal injuries or damage to the instrument. The PreON must only be operated by qualified personnel who have been appropriately trained.

Servicing of the PreON must only be performed by a PreOmics Field Service specialist.

WARNING

Risk of personal injury and material damage

[W4]

Do not use damaged rotor adapters. The rotor adapters can only be used once. High g forces exerted in the centrifuge can cause damage to used rotor adapters.

CAUTION

Damage to the instrument

[C3]

Empty the tip disposal container prior to use to prevent a tip jam in the waste drawer. Failure to empty the waste container may block the robotic arm that could cause run failure or instrument damage.

WARNING

Risk of personal injury and material damage

[W5]

To prevent plastic crash, load the tubes properly. After a plastic crash, sharp plastic particles, toxic and aggressive reagents could be inside the centrifuge. Be careful when handling items inside the centrifuge and thoroughly clean the centrifuge according to the maintenance procedure.

CAUTION

Damage to the instrument

[C4]

Only use the correct volume of liquids.

Exceeding the recommended volume of liquids may damage the centrifuge rotor or instrument.

WARNING

Risk of fire or explosion

[W6]

When using flammable liquids on the PreON (Acetonitrile, Isopropyl alcohol), handle such liquids carefully and in accordance with the required safety regulations. If liquid has been spilled, wipe it off and leave the PreON hood open to allow flammable vapors to disperse.

WARNING

Risk of explosion

[W7]

The PreON is intended for use with reagents and substances supplied with PreOmics kits or other than outlined in respective Information for use. Use of other reagents and substances may lead to fire or explosion.

CAUTION

Damage to the instrument

[C5]

Do not lean against the touchscreen when it is pulled out.

WARNING

Samples containing infectious agents

[W12]

Some samples used with this instrument may contain infectious agents. Handle such samples with the greatest of care and in accordance with the required safety regulations.

Always wear safety glasses, 2 pairs of gloves, and a lab coat.

The responsible body (e.g., laboratory manager) must take the necessary precautions to ensure that the surrounding workplace is safe, and that the instrument operators are suitably trained and not exposed to hazardous levels of infectious agents as defined in the applicable Safety Data Sheets (SDSs), OSHA⁶ or ACGIH⁷ documents.

Venting for fumes and disposal of wastes must be in accordance with all national, state, and local health and safety regulations and laws.

⁶ OSHA : Occupational Safety and Health Administration (United States of America).

⁷ ACGIH : American Conference of Government Industrial Hygienists (United States of America).

WARNING

Hazardous chemicals

[W13]

Some chemicals used with this instrument may be hazardous or may become hazardous after completion of the protocol run.

Always wear safety glasses, gloves, and a lab coat.

The responsible body (e.g., laboratory manager) must take the necessary precautions to ensure that the surrounding workplace is safe and that the instrument operators are not exposed to hazardous levels of toxic substances (chemical or biological) as defined in the applicable Safety Data Sheets (SDSs), OSHA or ACGIH[†] documents.

Venting for fumes and disposal of wastes must be in accordance with all national, state, and local health and safety regulations and laws.

WARNING

Vapors from volatile solvents

[W14]

Some of the solvents used with the instrument are volatile. Use the instrument only with adequate ventilation according to the applicable Safety Data Sheets (SDS), OSHA or ACGIH[†] documents.

WARNING

Moving parts

[W17]

Avoid contact with moving parts during operation of the PreON. Under no circumstances should you place your hands under the robotic arm when it is lowering. Do not attempt to move any tip racks or tubes whilst the instrument is operating.

WARNING

Hot surface

[W21]

The shaker can reach temperatures of up to 83°C (181.4°F). Avoid touching it when it is hot. Carefully remove the samples after a run.

**WARNING/
CAUTION**

Risk of personal injury and material damage

[W22]

Only perform maintenance that is specifically described in this user manual.

5.1 PreON software

The PreON is operated through a touchscreen, which guides you step-by-step through the correct loading of the work deck and selection of the protocol.


General functions of the PreON touchscreen are described below.

Button/Icon	Function
	Enables the user to scroll up through a list.
	Enables the user to scroll down through a list.
	The software automatically proceeds to the next screen.
	Returns to the previous screen.
	Returns to previous screen without saving any changes.
	Enables the user to change certain settings (e.g., to edit a user account).
	Enables the user to delete certain settings (e.g., to delete a user).
Text fields	Enables to edit text or value. A pop-up keyboard enables these changes.
Row in tables	Can be pressed to select the respective row. Either item will be selected, or row will be highlighted.
	Press to show additional information to respective item
	Press to show important information that must be followed during run setup for respective item.
	Press to show additional information to respective item
	Press to show important information that must be followed during run setup for respective item.
	Navigation back to the home screen
	Tools/Maintenance functions
	Configuration
	Log out from the instrument


Note: The instrument's touchscreen does not support swiping gestures.

5.1.1 Entering text and numbers


To enter text or numbers, touch the respective field. The corresponding on-screen keyboard will be displayed.

Keyboard for entering a password

Keyboard for entering a sample name

Keyboard for editing a protocol parameter

For protocol parameters, the value range is displayed. In the example shown in the screenshot above, values from 50–100 µl can be entered, but only in steps of 10 µl.

Buttons and icons in the on-screen keyboard are described below.

Button/Icon	Function
	Remove character to the left.
	Clear all from the field.
	To type the next letter in upper-case. After the letter has typed, the keyboard will show lower-case letters again.
	Switch to upper-case letters. Allows to type a number of upper-case letters. Press the symbol again to return to lower-case letters.
	Show special characters.
	Return to letters.

Button/Icon	Function
OK	Confirm and close.
Cancel	Discard and close.

If the entered value is not correct, the field border will change to red. It is not possible to proceed to the next screen in this case. Press the field again and correct the value according to the range displayed next to the field.

5.2 Switching on the instrument and logging in

1. Close the instrument hood.
2. Switch on the instrument by pressing the power switch to the inner position. You will hear a sound (if sound settings are enabled), and a startup screen appears. The instrument automatically performs initialization tests. If the centrifuge lid was closed, it will open.

After the initialization is complete, the **Login** screen appears.

PREOMICS

Please log in to start

User ID

Jane


Password

Log in

Login screen

3. Enter **User ID** and **Password** using the on-screen keyboard.
4. Press **Log in**.
5. The **Home** screen appears.

In case of a failed login, an exclamation icon (!) and information screen will appear. Touch the respective field to enter the **User ID** and **Password** again, ensuring that you enter the information correctly. The **User ID** is case sensitive.

Information screen from failed login, for example, due to entering incorrect password

5.3 Setting up a protocol run

Protocol setup starts from the **Home** screen.


PreOmics standard protocols are installed on the PreON upon delivery.

Important: Before starting any protocol, thoroughly read the relevant PreOmics kit handbook.

The touchscreen software will guide you through the protocol run setup and loading steps. The display screens vary depending on the protocol in use and may look different from the screens shown in this section.

1. Prepare the samples to be processed using the protocol indicated in the protocol selection screen. Refer to the kit handbook for more information, if required.

- To start the run setup, select from **Applications**, **Favorites** or **Scan kit**. To scan the 2-D bar code of the PreOmics kit, press **Scan Kit** and then use the handheld scanner, or just scan the bar code.

The software will automatically proceed to the next screen. If you scanned a kit bar code or if you selected a protocol from **Favorites**, the software may skip the **Kit**, **Material** and/or **Protocol Selection** screens. The software will skip the selection screen if the required information will be provided by kit bar code scanning or selection from **Favorites**.

- To enter information in the next screens, follow the instructions provided in the sections below. Depending on your selections, the number and order of the screens displayed on your instrument may vary.

Each section below contains a screenshot image. Follow the instructions in that section with the corresponding screen displayed on your instrument.

In general, press **Next** to proceed to the next screen, or press **Back** to return to the previous screen. **Next** will only be active if all required information has been entered in the current screen.

iST-NHS protocol

Step 1
Scan barcode chemical labeling reagent

PREOMICS

SetupToolsConfigurationLogout

Scan barcode chemical labeling reagent

Kit 1:

Clear

Cancel

BackNext

May 02, 2022, 13:21

User: Admin Admin

NHS

Step 1/17

Please scan the chemical label kit barcode. The kit barcode will be listed in the run report.

Step 2
Select sample type

PREOMICS

SetupToolsConfigurationLogout

Select sample type

☒ Pellet

☐ Liquid concentrated1 µl

☐ Liquid diluted11 µl

☐ Off-board lysis50 µl

Cancel

BackNext

May 02, 2022, 13:22

User: Admin Admin

NHS

Step 2/17

If sample type 'Pellet' is chosen PreON adds 50µl LYSE.

Step 3
Select digestion duration

PREOMICS

SetupToolsConfigurationLogout

Select digestion duration

☒ 1 hour

☐ 2 hours

☐ 3 hours

Cancel

BackNext

May 02, 2022, 13:22

User: Admin Admin

NHS

Step 3/17

Choose the duration of the digestion step. PreOmics recommends 1 hour for precipitated proteins, cells and body fluids. For tissue samples, 3 hours are recommended.

Step 4
Select number of samples

PREOMICS

SetupToolsConfigurationLogout

Select number of samples

NHS

Step4/17

☐ 4 samples

☒ 8 samples

☐ 12 samples

☐ 16 Samples

For less samples than indicated, use blank samples.

Cancel

BackNext

May 02, 2022, 13:23

User: Admin Admin

Step 5
Activate Low Input protocol

PREOMICS

SetupToolsConfigurationLogout

Low input protocol

NHS

Step5/17

☐ Activate Low Input protocol

Activate the low Protein input protocol to process 20 µg protein starting material or less.

The PreON will automatically adjust the volume of LYSE and DIGEST to 10 µl each, while keeping all other buffers in line with our protocol.

Cancel

BackNext

May 02, 2022, 13:24

User: Admin Admin

Step 6
Select chemical labeling plex degree

PREOMICS

SetupToolsConfigurationLogout

Select chemical labeling plex degree

NHS

Step6/17

☒ 2-plex

☐ 6-plex with 2 blanks

Choose the plex degree for the chemical labeling.

Cancel

BackNext

May 02, 2022, 13:27

User: Admin Admin

Step 7
Select amount of chemical
labeling reagent

PREOMICS

SetupToolsConfigurationLogout

Select amount of chemical labeling reagent

NHS
Step7/17

☐ 200 µg

☐ 600 µg

☒ 400 µg

☐ 800 µg

Select amount of chemical labeling reagent per channel. For example, if 50 µg of sample are used, 400 µg labeling reagent are required for a labeling ratio of 1:8.

Cancel

BackNext

May 02, 2022, 13:29User: Admin Admin

Step 8
Chemical labeling options

PREOMICS

SetupToolsConfigurationLogout

Chemical labeling options

NHS
Step8/17

Samples	Chemical label	Info
1, 2	126	
3, 4	126	
5, 6	127	
7, 8	127	

Edit the chemical label name if necessary. The name will be added to the sample ID and will be listed in the run report.

Cancel

BackNext

May 02, 2022, 13:31User: Admin Admin

Step 9
Select labeling time

PREOMICS

SetupToolsConfigurationLogout

Select labeling duration

NHS
Step9/17

Labeling time60min

The allowed labeling time is 60 - 120 minutes.

Cancel

BackNext

May 02, 2022, 13:33User: Admin Admin

Step 10
Select split or complete protocol

PREOMICS Setup Tools Configuration Logout

Select split or complete protocol NHS Step 10/17

☒ Split Protocol (recommended)
☐ Complete Protocol

Select 'Split Protocol' for best 'Chemical Labeling' performance. If this setting is chosen, the 'Chemical Labeling' reagents have to be loaded during the run. The PreON pauses after the digestion step.

Cancel Back Next

May 02, 2022, 13:35 User: Admin Admin

Step 11
Load ELUTE later

PREOMICS Setup Tools Configuration Logout

Load ELUTE later NHS Step 11/17

☒ Load ELUTE later (recommended)

If checked, the PreON will stop prior to the timepoint in the workflow where the ELUTE buffer is needed. Please insert the ELUTE buffer at this timepoint and then proceed as the PreON will indicate. For digestion times > 1h and in general for labelling experiments we recommend to check this option.

Cancel Back Next

May 02, 2022, 13:37 User: Admin Admin

Step 12
Sample drying

PREOMICS Setup Tools Configuration Logout

Sample drying NHS Step 12/17

☐ Enable Sample Drying

Enable drying of samples after peptide elution from cartridges. Duration is 4 hours for 4 and 8 samples, and 5 hours for 12 samples.

PreON will add defined LC-LOAD volume to sample after drying, we recommend a final concentration of 1 µg/µL.

Cancel Back Next

May 02, 2022, 13:38 User: Admin Admin

Step 13
Review Settings

PREOMICS

SetupToolsConfigurationLogout

Review

NHS
Step13/17

Sample typePellet

Digest duration60 min

Nr. of samples8

Low inputNo

Plex degree2-plex

AcetonitrilIn: 80 µl out: 40 µl

Labeling time60 min

PoolingNo

ProtocolSplit

Load ELUTE laterYes

Sample dryingNo

Review of your defined protocol, please check the conditions applied.

Press next to proceed or back to adjust the chosen settings.

Cancel

BackNext

May 02, 2022, 13:39

User: Admin Admin

Step 14
Insert 1

PREOMICS

SetupToolsConfigurationLogout

Insert 1

NHS
Step14/19

Put the insert Reagent 32 into slot "Insert 1".

Cancel

Move arm leftMove arm right

BackNext

May 02, 2022, 13:41

User: Admin Admin

Step 15
Insert 3

PREOMICS

SetupToolsConfigurationLogout

Load reagents

NHS
Step16/19

Pos.	Name	Volume	Info
B1	WASH 2	1925 µl	i
C1	WASH 1	1850 µl	i
D1	STOP	1050 µl	i
E1	Digest Tube (0.5 ml)	0 µl	i
E2	Digest Tube (0.5 ml)	0 µl	i
F1	RESUSPEND	675 µl	i
G1	LYSE	650 µl	i

Insert 1

A1A4

B1B4

C1C4

D1D4

E1E4

F1F4

G1G4

Cancel

Move arm leftMove arm right

BackNext

May 02, 2022, 13:43

User: Admin Admin

Step 18-2
Load the centrifuge rotor
adapter – Middle Position

PREOMICS

SetupToolsConfigurationLogout

Load rotor adapters (wash spin column)

NHS
Step18/19

Pos.	Labware	Lid position
1	PreOmics cartridge inside rotor cartridge adapter	-
2	empty	-
3	1.5 ml collection tube	L3

8 x Rotor adapter

L1

1

2

3

rotor inside

rotor outside

Change View

CancelMove arm leftMove arm rightBackNext

May 02, 2022, 13:49

User: Admin Admin

Step 18-3
Load the centrifuge rotor
adapter – Elute Position

PREOMICS

SetupToolsConfigurationLogout

Load rotor adapters (wash spin column)

NHS
Step18/19

Pos.	Labware	Lid position
1	PreOmics cartridge inside rotor cartridge adapter	-
2	empty	-
3	1.5 ml collection tube	L3

8 x Rotor adapter

L1

1

2

3

rotor inside

rotor outside

Change View

CancelMove arm leftMove arm rightBackNext

May 02, 2022, 13:52

User: Admin Admin

Step 18-4
Load the centrifuge rotor
adapter

PREOMICS

SetupToolsConfigurationLogout

Load rotor adapters (wash spin column)

NHS
Step18/19

Pos.	Labware	Lid position
1	PreOmics cartridge inside rotor cartridge adapter	-
2	empty	-
3	1.5 ml collection tube	L3

8 x Rotor adapter

L1

1

2

3

rotor inside

rotor outside

Change View

CancelMove arm leftMove arm rightBackNext

May 02, 2022, 13:53

User: Admin Admin

Step 19
Load shaker

PREOMICS

Setup

Tools

Configuration

Logout

Load shaker

NHS

Step 19/19

Pos	Sample ID	Tube Type	Volume
01	202205...01_126	1.5 ml	0 µl
02	202205...02_126	1.5 ml	0 µl
03	202205...03_126	1.5 ml	0 µl
04	202205...04_126	1.5 ml	0 µl
07	202205...05_127	1.5 ml	0 µl
08	202205...06_127	1.5 ml	0 µl
09	202205...07_127	1.5 ml	0 µl

13

01

02

03

04

05

06

07

08

09

10

11

12

16

Cancel

Move arm left

Move arm right

Back

Next

May 02, 2022, 13:54

User: Admin Admin

iST protocol

Step 1 Select sample type

PREOMICS Setup Tools Configuration Logout

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☒ Pellet
☐ Liquid concentrated 1 µl
☐ Liquid diluted 11 µl
☐ Off-board lysis 50 µl

If sample type 'Pellet' is chosen PreON adds 50µl LYSE.

Cancel Back Next

May 03, 2022, 09:35 User: Admin Admin

Step 2 Select digestion duration

PREOMICS Setup Tools Configuration Logout

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☒ 1 hour
☐ 2 hours
☐ 3 hours

Choose the duration of the digestion step. PreOmics recommends 1 hour for precipitated proteins, cells and body fluids. For tissue samples, 3 hours are recommended.

Cancel Back Next

May 03, 2022, 09:36 User: Admin Admin

Step 3 Select number of samples

PREOMICS Setup Tools Configuration Logout

Select number of samples iST Step 3/11

☐ 4 samples
☒ 8 samples
☐ 12 samples

For less samples than indicated, use blank samples.

Cancel Back Next

May 03, 2022, 09:36 User: Admin Admin

Step 4
Low input protocol

PREOMICS

SetupToolsConfigurationLogout

Low input protocol

IST Step 4/11

☐ Activate Low Input protocol

Activate the low Protein input protocol to process 20 µg protein starting material or less.

The PreON will automatically adjust the volume of LYSE and DIGEST to 10 µl each, while keeping all other buffers in line with our protocol.

Cancel

BackNext

May 03, 2022, 09:38

User: Admin Admin

Step 5
Load ELUTE later

PREOMICS

SetupToolsConfigurationLogout

Load ELUTE later

IST Step 5/11

☐ Load ELUTE later (not needed)

If checked, the PreON will stop prior to the timepoint in the workflow where the ELUTE buffer is needed. Please insert the ELUTE buffer at this timepoint and then proceed as the PreON will indicate. For digestion times > 1h and in general for labelling experiments we recommend to check this option.

Cancel

BackNext

May 03, 2022, 09:39

User: Admin Admin

Step 6
Sample drying

PREOMICS

SetupToolsConfigurationLogout

Sample drying

IST Step 6/11

☐ Enable Sample Drying

Enable drying of samples after peptide elution from cartridges. Duration is 4 hours for 4 and 8 samples, and 5 hours for 12 samples.

PreON will add defined LC-LOAD volume to sample after drying, we recommend a final concentration of 1 µg/ µL.

Cancel

BackNext

May 03, 2022, 09:41

User: Admin Admin

Step 7
Review

PREOMICS

>>> Setup

Tools

Configuration

Logout

Review

IST Step 7/11

Sample typePellet

Digest duration60 min

Nr. of samples8

Low inputNo

Load ELUTE laterNo

Sample dryingNo

Review of your defined protocol, please check the conditions applied.

Press next to proceed or back to adjust the chosen settings.

Cancel

Back

Next

May 03, 2022, 09:42

User: Admin Admin

Step 8
Insert 1

PREOMICS

>>> Setup


Tools

Configuration

Logout

Insert 1

IST Step 8/13

Put the insert Reagent 32 into slot "Insert 1".

Cancel

Move arm left

Move arm right

Back

Next

May 03, 2022, 09:43

User: Admin Admin

Step 9
Insert 3

PREOMICS

>>> Setup


Tools

Configuration

Logout

Insert 3

IST Step 9/13

Put the insert TMT-11 into slot "Insert 3".

Cancel

Move arm left

Move arm right

Back

Next

May 03, 2022, 09:44

User: Admin Admin

Step 10
Load reagents

PREOMICS

SetupToolsConfigurationLogout

Load reagents

Pos.	Name	Volume	Info
A1	ELUTE	1050 µl	i
A2	ELUTE	1050 µl	i
B1	WASH 2	1925 µl	i
C1	WASH 1	1850 µl	i
D1	STOP	1050 µl	i
E1	Digest Tube (0.5 ml)	0 µl	i
E2	Digest Tube (0.5 ml)	0 µl	i

iST Step 10/13

Insert 1

A1

B1

C1

D1

E1

F1

G1

A4

B4

C4

D4

E4

F4

G4

Cancel

Move arm left

Move arm right

Back

Next

May 03, 2022, 09:45

User: Admin Admin

Step 11
Load tip racks

PREOMICS

SetupToolsConfigurationLogout

Load tip racks

Pos.	Name	Amount
S2	Tip Rack, 200 µl	Min. 23

iST Step 11/13

Insert 2

S1

S2

Cancel

Move arm left

Move arm right

Back

Next

May 03, 2022, 09:46

User: Admin Admin

Step 12-1
Load the centrifuge rotor
adapter – Spin Column

PREOMICS

SetupToolsConfigurationLogout

Load rotor adapters (wash spin column)

Pos.	Labware	Id position
1	PreOmics cartridge inside rotor cartridge adapter	-
2	empty	-
3	1.5 ml collection tube	L3

8 x Rotor adapter

L1

L2

L3

rotor inside

rotor outside

iST Step 12/13

1

2

3

Change View

Cancel

Move arm left

Move arm right

Back

Next

May 03, 2022, 09:47

User: Admin Admin

Step 12-2
Load the centrifuge rotor
adapter – Middle Position

PREOMICS

SetupToolsConfigurationLogout

Load rotor adapters (wash spin column)

IST Step 12/13

Pos.	Labware	Lid position
1	PreOmics cartridge inside rotor cartridge adapter	-
2	empty	-
3	1.5 ml collection tube	L3

8 x Rotor adapter

rotor inside rotor outside

Change View

CancelMove arm leftMove arm rightBackNext

May 03, 2022, 09:48

User: Admin Admin

Step 12-3
Load the centrifuge rotor
adapter – Elute Position

PREOMICS

SetupToolsConfigurationLogout

Load rotor adapters (wash spin column)

IST Step 12/13

Pos.	Labware	Lid position
1	PreOmics cartridge inside rotor cartridge adapter	-
2	empty	-
3	1.5 ml collection tube	L3

8 x Rotor adapter

rotor inside rotor outside

Change View

CancelMove arm leftMove arm rightBackNext

May 03, 2022, 09:50

User: Admin Admin

Step 12-4
Load the centrifuge rotor
adapter

PREOMICS

SetupToolsConfigurationLogout

Load rotor adapters (wash spin column)

IST Step 12/13

Pos.	Labware	Lid position
1	PreOmics cartridge inside rotor cartridge adapter	-
2	empty	-
3	1.5 ml collection tube	L3

8 x Rotor adapter

rotor inside rotor outside

Change View

CancelMove arm leftMove arm rightBackNext

May 03, 2022, 09:51

User: Admin Admin

Step 13
Load shaker

PREOMICS

Setup

Tools

Configuration

Logout

Load shaker

IST
Step 13/13

13

01

02

03

04

05

06

07

08

09

10

11

12

16

Pos	Sample ID	Tube Type	Volume
01	202205...943-01	1.5 ml	0 µl
02	202205...943-02	1.5 ml	0 µl
03	202205...943-03	1.5 ml	0 µl
04	202205...943-04	1.5 ml	0 µl
07	202205...943-05	1.5 ml	0 µl
08	202205...943-06	1.5 ml	0 µl
09	202205...943-07	1.5 ml	0 µl

Cancel

Move arm left

Move arm right

Back

Next

May 03, 2022, 09:53

User: Admin Admin

iST-BCT protocol

Step 1
Select sample type

PREOMICS

SetupToolsConfigurationLogout

Select sample type

BCT Step 1/11

☒ Pellet

☐ Liquid concentrated 1 µl

☐ Off-board lysis 50 µl

If sample type 'Pellet' is chosen PreON adds 50µl LYSE.

Cancel

BackNext

May 04, 2022, 09:37

User: Admin Admin

Step 2
Select digestion duration

PREOMICS

SetupToolsConfigurationLogout

Select digestion duration

BCT Step 2/11

☒ 1 hour

☐ 2 hours

☐ 3 hours

Choose the duration of the digestion step. PreOmics recommends 1 hour for precipitated proteins, cells and body fluids. For tissue samples, 3 hours are recommended.

Cancel

BackNext

May 04, 2022, 09:39

User: Admin Admin

Step 3
Select number of samples

PREOMICS

SetupToolsConfigurationLogout

Select number of samples

BCT Step 3/11

☐ 4 samples

☒ 8 samples

☐ 12 samples

For less samples than indicated, use blank samples.

Cancel

BackNext

May 04, 2022, 09:39

User: Admin Admin

Step 4
Low input protocol

PREOMICS

SetupToolsConfigurationLogout

Low input protocol

BCT
Step4/11

☐ Activate Low Input protocol

Activate the low Protein input protocol to process 20 µg protein starting material or less.

The PreON will automatically adjust the volume of LYSE and DIGEST to 10 µl each, while keeping all other buffers in line with our protocol.

Cancel

BackNext

May 04, 2022, 09:40

User: Admin Admin

Step 5
Load ELUTE later

PREOMICS

SetupToolsConfigurationLogout

Load ELUTE later

BCT
Step5/11

☐ Load ELUTE later (not needed)

If checked, the PreON will stop prior to the timepoint in the workflow where the ELUTE buffer is needed. Please insert the ELUTE buffer at this timepoint and then proceed as the PreON will indicate. For digestion times > 1h and in general for labelling experiments we recommend to check this option.

Cancel

BackNext

May 04, 2022, 09:41

User: Admin Admin

Step 6
Sample drying

PREOMICS

SetupToolsConfigurationLogout

Sample drying

BCT
Step6/11

☐ Enable Sample Drying

Enable drying of samples after peptide elution from cartridges. Duration is 4 hours for 4 and 8 samples, and 5 hours for 12 samples.

PreON will add defined LC-LOAD volume to sample after drying, we recommend a final concentration of 1 µg/µL.

Cancel

BackNext

May 04, 2022, 09:42

User: Admin Admin

Step 7
Review

PREOMICS

>>> Setup

Tools

Configuration

Logout

Review

BCT
Step 7/11

Sample typePellet

Digest duration60 min

Nr. of samples8

Low inputNo

Load ELUTE laterNo

Sample dryingNo

Review of your defined protocol, please check the conditions applied.

Press next to proceed or back to adjust the chosen settings.

Cancel

Back

Next

May 04, 2022, 09:42

User: Admin Admin

Step 8
Insert 1

PREOMICS

>>> Setup


Tools

Configuration

Logout

Insert 1

BCT
Step 8/13

Put the insert Reagent 32 into slot "Insert 1".

Cancel

Move arm left

Move arm right

Back

Next

May 04, 2022, 09:43

User: Admin Admin

Step 9
Insert 3

PREOMICS

>>> Setup


Tools

Configuration

Logout

Insert 3

BCT
Step 9/13

Put the insert TMT-11 into slot "Insert 3".

Cancel

Move arm left

Move arm right

Back

Next

May 04, 2022, 09:44

User: Admin Admin

Step 10
Load reagents

PREOMICS

SetupToolsConfigurationLogout

Load reagents

Pos.	Name	Volume	Info
A1	ELUTE	1050 µl	i
A2	ELUTE	1050 µl	i
B1	WASH 2	1925 µl	i
C1	WASH 1	1850 µl	i
D1	STOP	1050 µl	i
E1	Digest Tube (0.5 ml)	0 µl	i
E2	Digest Tube (0.5 ml)	0 µl	i

BCT

Step10/13

Insert 1

A1

B1

C1

D1

E1

F1

G1

A4

B4

C4

D4

E4

F4

G4

Cancel

Move arm left

Move arm right

Back

Next

May 04, 2022, 09:44

User: Admin Admin

Step 11
Load tip racks

PREOMICS

SetupToolsConfigurationLogout

Load tip racks

Pos.	Name	Amount
S2	Tip Rack, 200 µl	Min. 23

BCT

Step11/13

Insert 2

S1

S2

Cancel

Move arm left

Move arm right

Back

Next

May 04, 2022, 09:45

User: Admin Admin

Step 12-1
Load the centrifuge rotor
adapter – Spin Column

PREOMICS

SetupToolsConfigurationLogout

Load rotor adapters (wash spin column)

Pos.	Labware	Id position
1	PreOmics cartridge inside rotor cartridge adapter	-
2	empty	-
3	1.5 ml collection tube	L3

8 x Rotor adapter

L1

L2

L3

rotor inside

rotor outside

BCT

Step12/13

1

2

3

Change View

Back

Next

Cancel

Move arm left

Move arm right

May 04, 2022, 09:46

User: Admin Admin

Step 12-2
Load the centrifuge rotor
adapter – Middle Position

PREOMICS

SetupToolsConfigurationLogout

Load rotor adapters (wash spin column)

BCT Step 12/13

Pos.	Labware	Lid position
1	PreOmics cartridge inside rotor cartridge adapter	-
2	empty	-
3	1.5 ml collection tube	L3

8 x Rotor adapter

L1

1

2

3

rotor inside

rotor outside

L2

L3

Change View

CancelMove arm leftMove arm rightBackNext

May 04, 2022, 09:46

User: Admin Admin

Step 12-3
Load the centrifuge rotor
adapter – Elute Position

PREOMICS

SetupToolsConfigurationLogout

Load rotor adapters (wash spin column)

BCT Step 12/13

Pos.	Labware	Lid position
1	PreOmics cartridge inside rotor cartridge adapter	-
2	empty	-
3	1.5 ml collection tube	L3

8 x Rotor adapter

L1

1

2

3

rotor inside

rotor outside

L2

L3

Change View

CancelMove arm leftMove arm rightBackNext

May 04, 2022, 09:47

User: Admin Admin

Step 12-4
Load the centrifuge rotor
adapter

PREOMICS

SetupToolsConfigurationLogout

Load rotor adapters (wash spin column)

BCT Step 12/13

Pos.	Labware	Lid position
1	PreOmics cartridge inside rotor cartridge adapter	-
2	empty	-
3	1.5 ml collection tube	L3

8 x Rotor adapter

L1

1

2

3

rotor inside

rotor outside

L2

L3

Change View

CancelMove arm leftMove arm rightBackNext

May 04, 2022, 09:48

User: Admin Admin

Step 13 Load shaker

Pos	Sample ID	Tube Type	Volume
01	202205...943-01	1.5 ml	0 µl
02	202205...943-02	1.5 ml	0 µl
03	202205...943-03	1.5 ml	0 µl
04	202205...943-04	1.5 ml	0 µl
07	202205...943-05	1.5 ml	0 µl
08	202205...943-06	1.5 ml	0 µl
09	202205...943-07	1.5 ml	0 µl

Buttons: Cancel, Move arm left, Move arm right, Back, Next

Status: May 04, 2022, 09:49 User: Admin Admin

WARNING

Risk of fire or explosion

[W6]

When using flammable liquids on the PreON (Acetonitrile, Isopropyl alcohol), handle such liquids carefully and in accordance with the required safety regulations. If liquid has been spilled, wipe it off and leave the PreON hood open to allow flammable vapors to disperse.

5.3.1 Loading reagents and tips

Important: When the **Load reagents** screen appears, the robotic arm will automatically move slowly – even when the instrument hood is open – so you can access all of the loading positions. Always stand clear of the instrument while the robotic arm is moving. Wait until the robotic arm has completed its movements before you start to load or unload reagents or tip racks. After you are finished loading and you proceed from this screen, the robotic arm will automatically move back to its original position (above tip rack position 3).

There is one tip rack that can be used on the PreON. A blue rack for 200 µl filter-tips. The instrument uses the notches on the filter-tip rack to identify the type. Only use tips designed for use with the PreON.

Important: Do not use damaged filter tips. Do not load damaged tip racks onto the worktable.

Waste drawer

Ensure that you empty the waste drawer containing used disposable labware before every run to avoid waste accumulation.

Press **Next** to proceed to loading the centrifuge or shaker rack, depending on the selected protocol.


Important: After loading and you proceed from the **Loading tips and enzymes** screen, the robotic arm will automatically move back to its original position (above tip rack position 3).

5.3.2 Loading the centrifuge

The **Loading centrifuge** screen guides you through setting up the required rotor adapters and centrifuge for the run. Make sure to read all essential and critical information before proceeding to the next step.

Rotor adapters can be placed into a rotor adapter holder, enabling convenient and easy preparation and loading of columns. Place the spin column-based cartridge adapters and spin column-based cartridges, tubes or samples into the appropriate positions in each rotor adapter as instructed by the software.

Place the lids into the correct lid position of the rotor adapter as indicated on your screen in the **Lid position** table column and the rotor adapter illustration. Make sure that the lids are pushed all the way down to the bottom of the slots on the sides of the rotor adapter. Incorrectly positioned lids can break off during centrifugation and cause the protocol run to crash.

Correctly loaded rotor adapter. **A** Rotor adapter is correctly loaded, and the 1.5 ml microcentrifuge tube lid is in the correct position; **B** Correctly loaded rotor adapter seen from the side.

Incorrectly loaded rotor adapter. **C** Rotor adapter is incorrectly loaded with a 1.5 ml microcentrifuge tube. The lid of the tube is not pushed all the way down to the bottom of the slot of the rotor adapter and could break off during centrifugation (compare with part A of the above figure); **D** Incorrectly loaded rotor adapter seen from the side (compare with part B of the above figure).

E

1.5 ml microcentrifuge
tube lid is in the wrong
slot of the rotor adapter

Rotor adapter is incorrectly loaded with a 1.5 ml microcentrifuge tube. The lid of the tube is positioned in the wrong slot of the rotor adapter. During column transfer, the lid of the spin column-based cartridge could crash onto the lid of the 1.5 ml microcentrifuge tube, causing the protocol run to crash.

WARNING

Risk of personal injury and material damage

[W5]

To prevent plastic crash, load the tubes properly. After a plastic crash, sharp plastic particles, toxic and aggressive reagents could be inside the centrifuge. Be careful when handling items inside the centrifuge and thoroughly clean the centrifuge according to the maintenance procedure.

WARNING

Risk of personal injury and material damage

[W4]

Do not use damaged rotor adapters. The rotor adapters can only be used once. High g forces exerted in the centrifuge can cause damage to used rotor adapters.

WARNING

Risk of personal injury and material damage

[W27]

Make sure that lids from spin column-based cartridges and 1.5 ml microcentrifuge tubes are in the correct position and pushed all the way down to the bottom of the slots on the sides of the rotor adapter. Incorrectly positioned lids can break off during centrifugation.

WARNING

Risk of personal injury and material damage

[W28]

Be sure the lid is completely removed from the spin column. Spin column-based cartridges with partially removed lids may not be removed properly from the rotor, causing the protocol run to crash.

Loading the centrifuge


Tubes with samples need to be loaded onto the shaker unit (refer to section 5.3.3 Loading the shaker), and the centrifuge must be prepared according the following description.

PREOMICS Setup Tools Configuration Logout

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Cancel Move arm left Move arm right Back Next

May 03, 2022, 09:51 User: Admin Admin

The number and tube positions of rotor adapters required for the protocol run will be shown in the table and illustration on your screen. The table shows how to load and position each rotor adapter. The **Pos.** column indicates the position in the rotor adapter, and the **Lid position** column indicates where to place the lid of a particular tube.

For each rotor adapter:

1. Load each tube/spin column-based cartridge adapter with spin column-based cartridge in its correct position as indicated in the table on the screen. Touch the table row to highlight the particular tube position in the illustration below the table.
2. Ensure that the tubes and spin column-based cartridge adapter with spin column-based cartridge are pushed firmly into the appropriate rotor adapter position.
3. Make sure that the lids are pushed all the way down to the bottom of the slots on the sides of the rotor adapter. Ensure that you place the lids in the correct lid positions.
4. Repeat steps 1–3 until all rotor adapters have been prepared.
5. Place the loaded rotor adapters into the centrifuge buckets as shown on the right side of the screen. For ease of use and high process safety, the rotor adapters only fit into the centrifuge buckets in one orientation.
6. Press **Next** to proceed with loading the samples into the shaker. Follow the instructions in the corresponding section below.

5.3.3 Loading the shaker

The **Load shaker** screen guides you through the loading of the shaker.

Sample tubes must be loaded on the shaker.

PREOMICS Setup Tools Configuration Logout

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Pos	Sample ID	Tube Type	Volume
01	202205...943-01	1.5 ml	0 µl
02	202205...943-02	1.5 ml	0 µl
03	202205...943-03	1.5 ml	0 µl
04	202205...943-04	1.5 ml	0 µl
07	202205...943-05	1.5 ml	0 µl
08	202205...943-06	1.5 ml	0 µl
09	202205...943-07	1.5 ml	0 µl

01 02 03 04 05 06 07 08 09 10 11 12 13 16

Cancel Move arm left Move arm right Back Next

May 03, 2022, 09:53 User: Admin Admin

Loading of shaker; lid positions filled with shaker rack plugs

In this step, the software displays the shaker positions, tubes and volume to be loaded in the table and in the schematic on the right.. Make sure to read all essential and critical information before proceeding to the next step.

1. If required, change the default **Sample IDs** in the respective fields using the on-screen keyboard. You can enter the value manually or scan the sample bar code using a handheld scanner. The sample ID is initially created using the YYYYMMDD-HHMM-no. format.


Note: Sample IDs are part of the run reports and can be part of log files and audit trail. They are not encrypted.

2. Prepare the correct tubes.
3. Load the tubes into the correct position of the shaker rack. Touch the table row to highlight the position in the schematic at the right.
4. Depending on the tube type, the tube lid must be placed into the slot next to the tube, as shown on the screen or indicated in the **Information** icon (i). Make sure that the lid/shaker rack plug is securely seated in the slot. Do not place a lid or shaker rack plug next to an empty shaker rack position.
5. Press **Next** to proceed to starting the run

Do not use 1.5 ml microcentrifuge tubes on the shaker. These microcentrifuge tubes can cause filter-tips to stick during sample transfer. Using this tube on the shaker can damage the pipetting system and can cause a centrifuge crash.


5.4 Starting the protocol run

A confirmation message will appear when the final step on the last setup screen has been completed.


1. Empty the liquid waste bottle
2. Empty the solid waste in the waste drawer
3. Close the hood to continue.
4. Press **Start** to begin the run. The screen will display the estimated run duration.


Important: Do not open the instrument hood during a run. In the event the hood is opened during a run, the run will stop. If the tip adapter has picked up a tip, it must be removed manually.

During the run, you can press the **Run Details** button to display the run steps. To return to the run view, press **Close**.

5.5 End of the protocol run

When the protocol run is completed, the elution position and content will be shown on the right side of the screen. Make sure that you follow the proper procedures for storing and handling the samples.

Perform regular maintenance as described in section 6.3 Regular maintenance procedure before starting the next run.

5.6 Stopping a protocol run

In case of emergency, a run can be stopped by pressing the **Stop Run** button.

If a run is stopped, perform the daily maintenance as described in section 6.4 Daily maintenance procedure, and ensure that no plastic parts are present in the centrifuge before starting the next run.

Note: If a protocol run is stopped, the run cannot be restarted; the samples must be processed manually (refer to section 7.3.1 Protocol interruption).

5.7 Saving run reports to the USB stick

Run reports are automatically saved on the instrument after each run. The number of run reports on the instrument is limited to 200. If this number is reached, the oldest run report will be replaced by the new run report.

To transfer run reports to the USB stick, proceed as follows:

1. Press the **Tools** icon.

2. Press the **Data exchange** tab. The number of available run reports is shown on the screen.

PREOMICS

Setup Tools Configuration Logout

Data exchange Run modules Maintenance Testing

Run reports available: 5

☐ All run reports
☒ Recent reports 5
☐ Delete reports after download

Download

Support Package Create and download

USB Connected

May 03, 2022, 09:59 User: Admin Admin

3. If not yet connected, connect the USB stick to one of the USB ports at the left of the touchscreen.
4. To save all available run reports to the USB stick, select **All run reports**. To save only the recent reports, select **Recent reports**. To enter the number of reports to be saved, touch in the **Recent reports** field.
5. If you would like to delete reports from the instrument after download, press **Delete reports after download**.
6. Press **Download** to save the reports to the USB stick.

A confirmation message that run reports are successfully saved to the USB stick will appear. The USB stick can be removed from the instrument.

Important: Do not remove the USB stick while the files are downloading. Wait until the download is completed.

5.8 Logging out

To log out from the instrument, follow these steps:

1. Press **Logout** at the top right of the screen.

2. To log out, confirm the message with **OK**. To stay logged in, press **Cancel**.
3. The **Login** screen appears.

6 Maintenance Procedures

**WARNING/
CAUTION**

Risk of personal injury and material damage

Only perform maintenance that is specifically described in this user manual.

[W22]

The following maintenance procedures must be carried out to ensure reliable operation of the PreON:

- Regular maintenance: after each protocol run
- Daily maintenance: after the last protocol run of the day
- Monthly maintenance: every month
- Periodic maintenance: when necessary; at least every 6 months

Optionally, these procedures may be performed to check and ensure the reliability of operation of the PreON:

- Tightness Test: ensures the tightness of the tip adapter (e.g., after O-Ring change)

The software provides step-by-step guidance under **Tools/Maintenance** for the maintenance procedures listed above, except for the regular maintenance.

Following these procedures ensures that the PreON is free from dust and liquid spills.

Select the cleaning agent according to the objective of the cleaning procedure, the sample material used and the downstream assay.

WARNING

Risk of fire or explosion

When using flammable liquids on the PreON (Acetonitrile, Isopropyl alcohol), handle such liquids carefully and in accordance with the required safety regulations. If liquid has been spilled, wipe it off and leave the PreON hood open to allow flammable vapors to disperse.

[W6]

Before using any cleaning or decontamination methods except those recommended by the manufacturer, users should check with the manufacturer that the proposed method will not damage the equipment.

6.1 Cleaning agents

The following disinfectants and detergents are recommended for cleaning the PreON.

Note: If you want to use disinfectants different from those recommended, ensure that their compositions are like those described below.

General cleaning of the PreON:

- Mild, alcohol-free Detergents (e.g., Mikrozyd® AF sensitive)
- 70% ethanol may be used for cleaning the metallic surfaces of the worktable only; it must not be used for cleaning the plastic inserts or the PreON hood

6.1.1 Disinfection

Mild, alcohol-free disinfectants, e.g. Mikrozyd® AF sensitive.

Ethanol-based disinfectants must not be used for disinfection of the plastic inserts and can only be used for disinfection of metallic surfaces, such as the worktable or inside the centrifuge: e.g., 25 g ethanol and 35 g 1-propanol per 100 g liquid or mikrozyd® AF liquid (Schülke & Mayr GmbH).

General instructions

- Do not use spray bottles to spray cleaning or disinfectant liquids onto surfaces of the PreON workstation. Spray bottles should be used only for items that have been removed from the workstation.
- If solvents or saline, acidic, or alkaline solutions are spilt on the PreON or if PreOmics buffers splash the instrument hood, wipe the spilt liquid away immediately.
- Follow manufacturer's safety instruction for handling cleaning agents.
- Follow manufacturer's instruction for soaking time and concentration of the cleaning agents. Immersing for longer than the recommended soak time can harm the instrument.
- Do not use alcohol or alcohol-based disinfectants to clean the PreON hood. Exposing the PreON hood to alcohol or alcohol-based disinfectants will cause surface cracking. Clean the PreON hood with distilled water only or a mild detergent.
- Take care that no liquid runs down the touchscreen. Liquid may be drawn through the dust protection sealing by capillary forces and cause malfunction of the display. To clean the touchscreen, moisten a soft lint-free cloth with water, ethanol, or a mild detergent and carefully wipe the display. Wipe dry with a paper towel.
- Before using any cleaning or decontamination methods except those recommended by the manufacturer, users should check with the manufacturer that the proposed method will not damage the equipment

CAUTION

Damage to the instrument

[C10]

Do not use bleach, solvents, or reagents containing acids, alkalis, or abrasives to clean PreON.

CAUTION

Damage to the instrument

[C11]

Do not use spray bottles containing alcohol or disinfectant to clean surfaces of the PreON. Spray bottles should be used only to clean items that have been removed from the worktables.

WARNING

Risk of fire

[W24]

Do not allow cleaning fluid or decontamination agents to come into contact with the electrical parts of the PreON.

WARNING

Risk of electric shock

[W10]

Do not open any panels on the PreON.

Risk of personal injury and material damage

Only perform maintenance that is specifically described in this user manual.

WARNING

Hazardous chemicals and infectious agents

[W16]

The waste may contain toxic material and must be disposed of properly. Refer to your local safety regulations for proper disposal procedures.

WARNING

Risk of personal injury and material damage

[W1]

Improper use of the PreON may cause personal injuries or damage to the instrument. The PreON must only be operated by qualified personnel who have been appropriately trained.

Servicing of the PreON must only be performed by a PreOmics Field Service specialist.

WARNING

Risk of explosion

[W23]

When cleaning the PreON with alcohol-based disinfectant, leave the PreON hood open to allow flammable vapors to disperse.

Only clean the PreON when worktable components have cooled down.

WARNING

Risk of fire or explosion

[W6]

When using flammable liquids on the PreON (Acetonitrile, Isopropyl alcohol), handle such liquids carefully and in accordance with the required safety regulations. If liquid has been spilled, wipe it off and leave the PreON hood open to allow flammable vapors to disperse.

WARNING

Vapors from volatile solvents

[W14]

Some of the solvents used with the instrument are volatile. Use the instrument only with adequate ventilation according to the applicable Material Safety Data Sheets (MSDSs) or OSHA,* ACGIH[†] or COSHH[‡] documents.

WARNING

Toxic fumes

[W14]

Do not use bleach to clean or disinfect the PreON. Bleach in contact with salts from the buffers can produce toxic fumes.

6.2 Servicing

Contact your local PreOmics Field Service Specialist or your local distributor for more information about flexible Service Support Agreements from PreOmics.

WARNING


Risk of personal injury and material damage

[W1]

Improper use of the PreON may cause personal injuries or damage to the instrument. The PreON must only be operated by qualified personnel who have been appropriately trained. Servicing of the PreON must only be performed by a PreOmics Field Service specialist.

6.3 Regular maintenance procedure

After running a protocol, perform the regular maintenance procedure described below.

Regular maintenance procedure:

1. Open the waste drawer and empty tips and columns (if necessary) into a suitable laboratory waste container.
2. Remove used disposable labware and unwanted samples and reagents from the worktable. Discard them according to your local safety regulations.

Note: If the robotic arm prevents you from reaching a position, do not move the robotic arm manually. Instead, proceed as follows:

Press **Move left** or **Move right**, as needed. The robotic arm will start to move. The hood can remain open during this movement.

Ensure that you stand clear of the instrument while the robotic arm is moving. Wait until the robotic arm has completed its movements.

3. Replace the lids of the reagent bottles and close tightly. Store the bottles according to the instructions in the relevant kit handbook.

You can now run another protocol or switch off the PreON.


6.4 Daily maintenance procedure

After running the last protocol of the day, perform the daily maintenance procedure. The software guides you through each step to be performed:

1. To start the daily maintenance, press the **Tools** icon.

2. Then press the **Daily** subtab under the **Maintenance** tab. The screen shows the **Last Executed** and the **Next Due** daily maintenance dates.

3. Press **Start**. Follow the instructions on the screen. Details are provided in the next steps below.
The robotic arm will automatically move slowly to the left – even if the instrument hood is open – to provide access to the loading positions. Always stand clear of the instrument while the robotic arm is moving. Wait until the robotic arm has completed its movements before you start to unload.
4. Remove used disposable labware, adapters and unwanted samples and reagents from the worktable. If required, discard them according to your local safety regulations.
5. Close the buffer bottles tightly and store according to the instructions in the relevant kit handbook.
6. Press **Done** to confirm that the steps have been completed.
7. Empty the waste drawer and check that the inlay is clean. If necessary, clean the inlay of the waste drawer with alcohol-based disinfection wipes, or by soaking using one of the cleaning agents listed above, and then rinse with distilled water.
8. Wipe and clean the worktable with alcohol-based disinfection wipes. Incubate as appropriate, rinse thoroughly with distilled water and wipe dry with paper towels.
Note: Do not use alcohol or alcohol-based disinfectants to clean the hood.
9. Press **Done** only when the steps listed above have been successfully completed. The date of the last performed daily maintenance is updated automatically.
The robotic arm will automatically move back to its original position (above tip rack position 3).


6.5 Monthly maintenance procedure

Perform the daily maintenance procedure before you perform the monthly maintenance procedure. Select the appropriate cleaning agent according to the sample material and downstream assay (refer to section 6.1 Cleaning agents).

1. To start the monthly maintenance, press **the Tools icon**.

2. Then press the **Monthly** subtab under the **Maintenance** tab. The screen shows the **Last Executed** and the **Next Due** monthly maintenance dates.

3. Close the hood.
4. See Cleaning procedure section 6.1 Cleaning agents
5. Incubate the shaker adapter (grey), shaker tray (metal adapter), buffer bottle rack and waste drawer in liner by soaking as appropriate. Rinse thoroughly with distilled water and wipe dry with paper towels.
6. Press **Done** only when the steps listed above have been successfully completed. The date of the last performed monthly maintenance is updated automatically.
Important: Inspect the waste drawer during maintenance. Contact PreOmics Technical Services if any broken parts are observed.
7. Transfer the run reports from the instrument to the USB stick and remove the run reports from the instrument. Only the last 200 run reports are saved on the instrument. For details refer to section 5.7 Saving run reports.

6.6 Periodic maintenance procedure

The periodic maintenance consists of cleaning the robotic arm modules and centrifuge.


Select the appropriate cleaning agent according to the sample material and downstream assay (refer to section 6.1 Cleaning agents).

6.6.1 Cleaning the robotic arm modules

Cleaning of the robotic arm modules must be performed periodically or could be performed if required. For example, the robotic arm modules must be cleaned if liquids were spilled due to crash.

Note: Perform the monthly maintenance procedure before you perform the robotic arm cleaning procedure.

1. To start cleaning the robotic arm modules, press the **Tools** icon and press the **Robotic Arm** subtab under the **Maintenance** tab. The screen shows the **Last Executed** maintenance date of the robotic arm modules.

2. Press **Start** to begin the cleaning of robotic arm modules. Follow the instructions on the screen. Details are provided in the next steps below.
3. Make sure that used Labware, adapters and reagents are removed from the worktable. Close the hood.
4. Press **Next** to move to cleaning position.
5. Remove the waste drawer and open hood.
6. Open the waster drawer. Moisten a soft lint-free cloth with water and carefully clean the optical sensor, tip adapter, gripper unit, rotor adapter stabilization rod and the spin column-based cartridge lid holder. Wipe these items dry as indicated on the touchscreen of the instrument.
7. Close the hood and press **Done** to finish cleaning of robotic arm. The date of the last performed cleaning of robotic arm is updated automatically.

6.6.2 Cleaning the centrifuge

Cleaning of the centrifuge arm must be performed periodically or could be performed if required. For example, the centrifuge must be cleaned in case of plastic crash or spillage of liquids due to crash.

WARNING


Risk of personal injury and material damage

[W5]

To prevent plastic crash, load the tubes properly. After a plastic crash, sharp plastic particles, toxic and aggressive reagents could be inside the centrifuge. Be careful when handling items inside the centrifuge and thoroughly clean the centrifuge according to the maintenance procedure

Note: Perform the monthly maintenance procedure before you perform the cleaning of centrifuge procedure.

1. To start cleaning the centrifuge, press the **Tools** icon and press the **Centrifuge** subtab under the **Maintenance** tab. The screen shows the **Last Executed** centrifuge maintenance date.

2. Press **Start** to begin the centrifuge cleaning procedure. Follow the instructions on the screen. Details are provided in the next steps below.
3. The centrifuge lid must be open to allow access to the inside of the centrifuge. The lid should be opened only after the centrifuge has come to a complete stop. If the lid does not open automatically, close the hood and press the **Open Centrifuge Lid** button.
4. Perform cleaning as outlined in the following sections (below): Cleaning the rotor and buckets, Cleaning the centrifuge, Maintenance of the rotor nut, and Installing the centrifuge rotor and buckets.
5. Turn on the instrument. Press the **Tools** icon and press the **Centrifuge** subtab under the **Maintenance** tab.
6. When cleaning is completed, press **Start** again then press **Done** to confirm cleaning. The date of the last performed cleaning of centrifuge is updated automatically.

Cleaning the rotor and buckets

1. Switch off the PreON.
2. Remove all disposable rotor adapters, including tubes and spin column-based cartridges, from the buckets.
3. Remove the buckets from the rotor. Undo the rotor nut on top of the rotor using the rotor key, and carefully lift the rotor off the rotor shaft.

Rotor key

4. Submerge the rotor, buckets and rotor nut in cleaning agent. Incubate as appropriate.
5. Rinse thoroughly with distilled water. Use a brush (e.g., a toothbrush or tube brush) to clean any parts that are difficult to access, such as the bucket mount and the rotor head. Wipe surfaces dry with a soft lint-free cloth. If available, dry the buckets and rotor with pressurized air.

Brushing a bucket

Brushing the rotor

Important: Make sure the paper towels and brush used are lint-free.

Important: Make sure that all residual salt is removed.

Important: Make sure to remove all traces of cleaning agent from the centrifuge buckets. Residual agent can cause the buckets to jam.

6. Carefully check the rotor for damage. If the rotor is damaged or shows signs of wear or corrosion, do not use the rotor. Contact PreOmics Technical Services.
7. Apply a few drops of mineral oil (Anti-Corrosion Oil (rotor)) on a soft, lint-free cloth, and wipe the bucket mount and rotor claw. A thin, invisible oil film should cover the bucket mount and rotor claw, but no droplets or smear should be apparent.
8. Apply oil to the rotor claw and to the bucket mount.

Important: Before applying oil to the rotor buckets on the rotor, make sure that the rotor and all buckets are completely dry.

Rotor head

Bucket mounts

Cleaning the centrifuge

1. Moisten a soft lint-free cloth with cleaning agent and clean the inside of the centrifuge and the centrifuge gasket. Incubate as appropriate.
2. Clean the inside of the centrifuge and the gasket with distilled water and wipe dry with lint-free paper towels. If available, use a vacuum cleaner.

Important: Make sure the gaskets remain in the proper positions.

3. Clean the centrifuge lid with a soft lint-free cloth moistened with cleaning agent. Incubate as appropriate, clean with water and wipe dry with paper towels.

4. Check the centrifuge gasket for damage. If the gasket is damaged or shows signs of wear, contact PreOmics Technical Services.

Maintenance of the rotor nut


Note: Always perform the cleaning procedure after disassembly of the rotor and at least twice a year.

After cleaning the rotor thread, apply a few drops of mineral oil (Anti-Corrosion Oil (rotor)) on a lint-free cloth, and wipe the thread. A thin, invisible oil film should cover the rotor thread but no droplets or smear should appear.

Rotor
thread

Rotor thread.

Inner
thread

Inner thread of the rotor nut.

After cleaning the inner thread of the rotor nut, wipe the thread using Anti-Corrosion Oil as described above.

Installing the centrifuge rotor and buckets

1. Mount the rotor.
2. The rotor can be mounted in only one orientation. The pin on the rotor shaft fits into a notch on the underside of the rotor directly underneath rotor position 1. Line up position 1 of the rotor with the pin on the rotor shaft and carefully lower the rotor onto the shaft.
3. Install the rotor nut on top of the rotor and tighten securely using the rotor key supplied with the PreON. Make sure that the rotor is securely seated.

Rotor key.

Rotor nut.

If the rotor nut is not tightened properly, it can become loose during operation of the centrifuge and can cause serious damage to the instrument. Such damage is not covered by the warranty.

WARNING

Risk of personal injury and material damage

[W25]

To prevent the rotor nuts from loosening during operation of the centrifuge, securely tightened the nuts using the rotor key supplied with PreON.

4. Replace the rotor buckets.

When replacing the rotor buckets, the side of the rotor bucket that must face toward the rotor shaft is marked with a grey line. Hold the bucket at an angle with the grey line facing the center of the rotor and hang the bucket on the rotor. Check that all buckets are properly suspended and can swing freely.

Important: All centrifuge buckets must be mounted before starting a run.

Before starting next protocol run, the centrifuge must be operated independently before starting further runs to check if residual plastic parts are still in the centrifuge. Follow the instructions below:

1. Login as User
2. To start a centrifuge run, press the Tools icon and then the Run Modules tab.
3. In the Set speed and Set duration fields, set the speed to 10,000 g and the duration to 1 min (1:0 m:s), respectively.
4. Press Start to begin the centrifuge run.
5. Carefully listen to the sound during centrifugation. See below for more details regarding the sound.

Unusual sound during centrifugation

If any grinding, rattling or crunching sounds are heard during the centrifugation, there are still loose plastic particles inside the centrifuge. Repeat the cleaning procedure as described above.

Note: It may be necessary to repeat the procedure several times to remove all plastic particles.

No unusual sound during centrifugation

If no unusual sound from loose plastic particles can be heard during centrifugation, the next protocol run can be started.


Note: The Open centrifuge lid and Close centrifuge lid buttons are not needed to start a centrifuge run, as the lid will close automatically. Instead, they are needed in case you need to prepare the PreON for shipment or during troubleshooting.

6.6.3 Tightness test

To ensure that the tightness of the tip adapter is sufficient for accurate pipetting, the tightness test of the tip adapter must be performed. This test must also be performed after replacing a tip adapter O-Ring to verify if replacement is successful.

Note: Perform the daily maintenance and cleaning of the robotic arm procedure before you perform the tightness test.

- To start the tightness test, press the **Tools** icon and press the **Tightness** subtab under the **Maintenance** tab. The screen shows the **Last Executed** tightness test date.

- Press **Start** to begin with the tightness test procedure. Follow the instructions on the screen. Details are provided in the next steps below.
- Open the hood and load a 200 µl tip rack with at least one 200 µl tip into tip rack position S3.
- Place an empty 1.5 ml tube in position 1 of the shaker.
- Place a 2 ml screw-cap microtube filled with 2 ml 96–100% ethanol in position H1.
- Close the hood and press **Next** to start tightness test.
After the load check, the robotic arm will pick up a tip, aspirate ethanol and move to the tube. The tip will remain in place above the tube for 2 minutes. The tip will be discarded into the waste afterwards.
- Wait until the test has been completed and then press **Next**.
- Remove the tube and visually check if liquid is present:
If no liquid is present, press **Yes** to record that the test passed.
If liquid is present, press **No** to record that the test failed.
- In case the test failed, repeat the test. If test fails again, it is recommended to replace the O-Ring first (using the O-Ring tab) or contact PreOmics Technical Services.
- Remove all remaining consumables and store accordingly
- Press **Done** to finalize tightness test procedure. The date of the last performed tightness test is updated automatically.

6.6.4 O-Ring exchange

O-Ring replacement must be performed if the tightness test failed or if the following issues are observed:

- Uneven volume transfers
- Dripping on the worktable

The replacement procedure requires the O-Ring change tool and an O-Ring. See Appendix B – PreON Accessories for ordering details.


O-Ring tool with prepared new O-Ring

The O-Ring replacement is semi-automatically and includes movement of robotic arm.

Note: Perform the daily maintenance and cleaning of the robotic arm procedure before you replace the O-Ring.

1. To start the O-Ring replacement, press the **Tools** icon and press the **O-Ring** subtab under the **Maintenance** tab. The screen shows the last O-Ring replacement date.


2. Close the hood and press **Start** to begin the O-Ring replacement procedure. Follow the instructions on the screen. Details are provided in the next steps below.
3. To prepare the O-Ring tool, perform the following steps:
 - 3a. Slide the new O-Ring over the small end of the peg.


- 3b. Push the grey lever until you reach the black stop and insert the small end of the peg into the hole.
- 3c. Press the peg down using the back end of the tweezers until the O-Ring sits (in the middle) on the larger end of the peg.

- 3d. Open the grey lever and insert the peg with small end first into hole as shown.

- 4. Press **Next** to begin loading of the O-Ring tool into the PreON.
- 5. Load the O-Ring tool by opening the grey lever into tip rack position 1 (nearest to user).
- 6. Close the hood and press **Next** to begin cutting of the O-Ring.
- 7. To cut and remove the O-Ring, perform the following steps:
 - 7a. To cut the O-Ring, open the hood and rotate the grey lever counter-clockwise until you reach the black stop.

- 7b. Open the grey lever and remove the O-Ring (by using the tweezers) from pipetting channel.
Note: If required, repeat cutting process until O-Ring is cut completely and can be removed.

8. Close the hood and press **Next** to pick up the prepared new O-Ring.
9. Open the hood and visually check if new O-Ring sits firmly on tip adapter.

- Note:** If the O-Ring was not successfully picked up, complete the O-Ring replacement procedure and restart.
10. Close the hood press **Next** to remove the O-Ring change tool.
11. Open the hood and remove the O-Ring change tool.

12. Wipe and clean the O-Ring change tool with alcohol-based disinfection wipes. Incubate as appropriate, rinse thoroughly with distilled water and wipe dry with paper towels.
13. Press **Done** to complete the O-Ring replacement. The date of the last performed O-Ring replacement is updated automatically.

WARNING

Risk of personal injury and material damage

[W1]

Improper use of the PreON may cause personal injuries or damage to the instrument. The PreON must only be operated by qualified personnel who have been appropriately trained.

Servicing of the PreON must only be performed by a PreOmics Field Service specialist.

6.7 Decontaminating the PreON

If the PreON is contaminated with infectious material, it should be decontaminated. If hazardous material is spilt on or inside the PreON, the user has responsibility for carrying out appropriate decontamination.

The PreON has to also be decontaminated before handling by a service engineer or shipping. In this case, a decontamination certificate must be completed to confirm that the decontamination procedure has been carried out.

To decontaminate the PreON, follow the daily, monthly and periodic maintenance procedure in sections 6.4–6.6.4, using the recommended disinfection agents.

7 Troubleshooting

7.1 General information

This section provides information about what to do if an error occurs while using the PreON.

7.2 Contacting PreOmics Technical Services

Whenever encountering a PreON error, be sure to have the following information at hand:

- Protocol name and version (found in the report file)
- Software version (see section 4.5.6).
- Sample input material
- Detailed description of the error situation

This information will help you and your PreOmics Technical Service Specialist to deal most efficiently with the issue. Please reach out to your local PreOmics contact for assistance. If unsure, get in touch at preomics.com/contact

In some cases, updates may be available for addressing specific problems.

7.2.1 Creating a support package

1. On the main screen, press the **Tools** button.

2. Press the **Data Exchange** tab.
3. Connect the USB stick to one of the 2 USB ports next to the touchscreen.

4. Press **Create and Download**. The support package will be created and saved on the USB stick in the folder **Support_Package**. The support package is a zip file that can be sent to PreOmics Technical Services.

7.3 Operation

	Comments and suggestions
Centrifuge	
Bucket does not swing back into place	Clean the centrifuge and rotors as described in section 6.
Imbalance detected	Make sure the rotor is symmetrically loaded according to instructions on run setup screens. Remove the rotor and check the centrifuge chamber for loose plasticware. Power OFF the PreON, wait for a few minutes, and power it ON again. If the error persists, contact PreOmics Technical Services.
Imbalance detected; loud noise heard during centrifugation	Ensure that loose parts are cleared from the worktable prior to centrifuge operation to avoid loose parts jamming or damaging the centrifuge.
Shaker	
Incorrect repositioning of shaker	The shaker should re-position itself towards the right side once shaking is completed. Remove any obstructions that prevent the shaker from returning to the correct position.
Robotic arm	
Robotic arm does not return to set position	Ensure that the instrument is placed on a stable, flat and level surface as described in section 4.1.1. In other cases, contact PreOmics Technical Services.
Pipettor	

	Comments and suggestions
Pipet tips not picked up by automatic pipettor	Make sure that the tip rack is not damaged and is correctly positioned on the worktable.
Pipet tips not disposed of correctly	Empty the waste drawer and ensure it is not broken. Check that tip disposal slot is not damaged or obstructed. Perform regular maintenance, as described in section 6.3
Droplets observed on worktable	The pipettor is dripping liquid. Check that the reagent bottles contain the correct buffers and are correctly placed in the reagent bottle rack. Be sure to use the right plasticware. Check the volumes in the sample tubes and tubes of accessory buffer(s), if applicable. Do not exceed the recommended amount of starting material to avoid blocking disposable filter-tips. If tip racks have been refilled, make sure that the correct tips were used. Check tightness of pipettor in maintenance section as described in section 6.6.3. If leakage could be detected change O-Ring as described in section 6.6.4. If issue persists, contact PreOmics Technical Services.
Mechanical	
Frame of instrument is distorted (e.g., uneven, unstable or not level)	Ensure that the instrument is placed on a stable, flat and level surface as described in section 4.1.1.
Hood sensor error: instrument will not function	Ensure the hood is properly closed. The instrument will not function if the hood is open.
Broken instrument hood	Ensure that only the cleaning products as described in section 6 are used on the hood.
Waste drawer jams but can still be inserted	Empty the waste drawer. Perform daily maintenance, as described in section 6.4.
Incorrectly inserted waste drawer	Handle the waste drawer with both hands when inserting or removing the drawer.
Pipet tips not disposed of correctly	Make sure the top of the tip disposal slot (refer to section 3.2) is not broken.
Scratches appear on the instrument	Always use the cleaning products as described in section 6. Do not use bleach or ethanol, as they can damage the surface of the instrument.
Electronic	
Display does not turn on	Do not touch the display with excessive force or use corrosive chemicals to clean the display surface. Contact PreOmics Technical Services for repair.
Error when copying files to USB	Power OFF the PreON, wait for a few minutes, and power it ON again. Save the file(s) to the USB stick again. Check the USB stick on a PC to ensure it is functional. If the error persists, contact PreOmics Technical Services.
USB device not detected	Power OFF the PreON, wait for a few minutes, and power it ON again. Insert the USB stick into the USB port. Check the USB stick on a PC to ensure it is functional. If the error persists, contact PreOmics Technical Services.

	Comments and suggestions
Login screen not visible when launching instrument	If the touchscreen does not display the login screen, but instead a software update message is shown, power OFF the PreON, wait for a few minutes. Ensure that the USB stick is not inserted in the USB port. Power ON the PreON again. The login screen should be visible. If the error persists, contact PreOmics Technical Services.
Error displayed when inserting the USB stick into a Windows PC	Ignore the message. In most cases, no scan is needed; use the USB stick as usual. Do not reformat the USB on the Windows PC. This will lead to complete data loss on the USB stick, and it can no longer be used with the PreON.

7.3.1 Protocol interruption

If an error occurs during a protocol run, it is possible to continue sample preparation manually. The error code, the description and the step at which the protocol stopped are displayed in the touchscreen.

To continue sample processing:

1. Note the step at which the protocol stopped. This is displayed in the touchscreen.
2. Remove the samples and reagents from the PreON.
3. Refer to the appropriate protocol in the relevant kit handbook and continue sample processing manually.

7.3.2 Centrifuge

Opening the centrifuge lid in the event of a breakdown

In case of power failure, the centrifuge lid can be manually opened so that the samples can be removed. To open the centrifuge lid, follow the instructions below.

WARNING

Moving parts

[W19]

In case of breakdown caused by power failure, remove the power cord and wait 10 minutes before attempting to manually open the centrifuge lid.

WARNING


Risk of personal injury and material damage

[W20]

Raise the centrifuge lid carefully. The lid is heavy and may cause injury if it falls.

1. Power OFF the PreON.
2. Unplug the power cord from the power outlet. Wait 10 minutes for the rotor to stop.
3. Open the instrument hood.
4. Carefully move the robotic arm to the right side of the worktable, furthest away from the centrifuge lid.

5. Remove the screw protection on top of the centrifuge lid. Using the rotor key, turn the screw counter-clockwise.

Turning the screw in the centrifuge lid

6. Remove the waste drawer. The centrifuge release cord will be visible on the left side of the waste drawer compartment.


Removed waste drawer

Centrifuge release cord

7. Pull the cord firmly to release the lid from the lock.
8. Manually raise the centrifuge lid.
9. Hold the raised lid and remove the samples and rotor adapters from the rotor.

Removing rotor adapters

Contact PreOmics Technical Services for instructions on how to reset the lid.

Liquid spills in the centrifuge

The rotor adapter is designed for use with PreOmics automated protocols. Do not fill the rotor adapters with liquid.

Liquid spills may occur if PreOmics spin column-based cartridges become blocked due to sample overloading. Do not exceed the maximum amount of starting material.

Incorrect installation of the centrifuge buckets may also cause rotor adapters to leak. Check that the buckets are installed properly and can swing freely.

If there is a liquid spill in the centrifuge, clean according to the instructions in section 6.

7.3.3 Reagent volume detection and ultrasonic pipe

The instrument does not start a load check if the ultrasonic pipe (black cap) of the ultrasonic sensor is missing. Check if cap has been installed before starting a load check.

Black beam columnator (see red circle) of the ultrasonic sensor

7.3.4 Touchscreen

Every time the user presses a button on the touchscreen, a small red sign is displayed at the place where the touchscreen recognizes the contact. If the point of touch and the recognized contact are at different positions, a re-calibration of the touchscreen can be performed. The calibration function can be reached during the instrument's startup procedure.

It is recommended to use a touch pen or an unused tip for optimal calibration results. In case you use a tip, discard the tip after calibration.

To re-calibrate the touchscreen:

1. Turn OFF the PreON.
2. Wait a few minutes and then turn ON the instrument again.


3. In the second screen press the PreOmics logo.

Note: If you do not press the logo, the instrument will continue initialization.

PREOMICS


PreON

4. Press **Calibrate** touchscreen.

5. An information message is displayed. The message will close automatically after 10 seconds.

6. Plus signs (+) will be shown at different positions on the screen. For each of these, press the center of the + sign. After one position is touched, the next position will be shown. The graphic below shows the positions and the order in which the + signs will appear.

7. The calibration is finished after all five positions have been touched. Calibration results are saved automatically. The previous screen is shown again.
8. Press **Start**. The instrument continues initialization using the new calibration settings.
- To cancel the calibration process, turn OFF the PreON.

8 Automated Status Supervision

The PreON has an internal web server which is accessible on port 80 with the HTTP protocol. The server delivers two different pages:


A HTML website with the status of the instrument.

/status.ini file with the status of the instrument in an ini-file format.

The pages are new generated on each call and reflect therefore always the actual state of the instrument.

8.1 HTML Website

The web site is quite simple:

8.2 /status.ini file

The purpose of the status.ini file is automated tracking of the instrument state. The file can be downloaded for example with `wget`:

```
wget HOSTNAME/status.ini
```

8.2.1 Format

Section	Name	Value
Preon	Serialnumber	The serial number of the instrument as a number.
PreonStatus	IsRunning	String "Yes" or "No".
	RunID	The run id as a number.
	Protocol	The name of the actual (or last) protocol (string).
	NumberOfSamples	The number of samples of the actual (or last) protocol (number).
	RunStatus	When running "Running". If the protocol ended without an error "Finished without error" or in case of an error "Finished with error (xx)".
	Duration	The duration of the actual (or last) protocol in [s].
	Remaining	The remaining time of the actual protocol in [s].
	UserIntervention	Empty string or in case of a required user intervention a describing text.
	Timestamp	Timestamp in ISO 8601 format.

Example:

```
[Preon]
Serialnumber=12
[PreonStatus]
IsRunning=No
RunID=0
Protocol=
NumberOfSamples=0
RunStatus=Finished without error
Duration=46
Remaining=17
UserIntervention=
Timestamp=2022-03-10T17:32:22
```

9 Glossary

Term	Description
Centrifuge	A component of the PreON that accommodates a rotor with 12 swing-out buckets. Each bucket holds a disposable rotor adapter.
Hood	The main door at the front of the PreON. When open, it provides complete access to the worktable.
Error code	A 3- or 4-digit number that indicates an error of the PreON.
Filter-tip	An item of labware that is picked up by the tip adapter during operation of the PreON. Liquid is aspirated into and dispensed from a filter-tip.
Initialization	An operation performed automatically when the PreON is switched on and if required before each protocol run to check the operation of the PreON.
Pipetting system/ Pipettor unit	The component of the PreON that aspirates and dispenses liquid. The pipetting system moves up and down above the worktable and contains a syringe pump that is connected to a tip adapter.
Power switch	A button located at the front of the PreON in the bottom-right corner. It allows the user to switch the PreON on and off; inner position is ON and outer position is OFF.
Protocol	A set of instructions for the PreON that allows the instrument to automate protein purification procedure. Protocols are run using the touchscreen.
Robotic gripper	A component of the PreON robotic arm that moves spin column-based cartridges during sample processing.
Rotor adapter	A disposable plastic adapter that fits into a centrifuge bucket and holds a PreOmics spin column-based cartridge and microcentrifuge tube during sample processing.
Tip adapter	A metal probe installed on the pipettor head. During operation of the PreON, the tip adapter picks up filter-tips from the worktable.
Tip disposal slots	Slots in the PreON worktable through which used filter-tips are discarded into the waste drawer.
Tip rack	A plastic rack that accommodates filter-tips on the worktable.
Touchscreen	The user interface that allows the user to operate the PreON.

Term	Description
Waste drawer	A drawer that collects used filter-tips and disposable columns.
Worktable	The surface of the PreON where samples, reagents, and filter-tips are loaded.

10 Appendix A – Technical data

10.1 Operating conditions

Power	100–240 V AC, 50/60 Hz, 650 VA Mains supply voltage fluctuations are not to exceed 10% of nominal supply voltages. Note: The apparent power can exceed 650 VA for up to 2 seconds during the centrifuge acceleration and can reach an approximate value of 1000 VA.
Fuse	2x T8A L 250V
Overvoltage category	II
Air temperature	18 to 28°C (64.4 to 82.4°F)
Relative humidity	15–75% (noncondensing)
Altitude	Up to 2000 m (6500 ft.)
Place of operation	For indoor use only
Pollution level	2
Environmental class	3K2 (IEC 60721-3-3)

10.2 Transport conditions

Air temperature	–25°C to 60°C (–13°F to 140°F) in manufacturer's package
Relative humidity	Max. 75% (noncondensing)
Environmental class	2K2 & 2M2 (IEC 60721-3-2)

10.3 Storage conditions

Air temperature	15°C to 30°C (59°F to 86°F) in manufacturer's package
Relative humidity	Max. 75% (noncondensing)
Environmental class	1K2 (IEC 60721-3-1)

10.4 Mechanical data and hardware features

Dimensions (hood closed)	Width: 65 cm (25.6 in.) Height: 58 cm (22.8 in.) Depth: 62 cm (24.4 in.)
Dimensions (hood open)	Width: 65 cm (25.6 in.) Height: 86 cm (34 in.) Depth: 62 cm (24.4 in.)

Mass	PreON: 71.5 kg (157.6 lb.) Accessories: 3 kg (6.6 lb.)
Centrifuge	10,640 rpm maximum 12,000 x g maximum Swing-out rotor, maximum 45° 12 rotor positions
Shaker	Speed 100–2000 RPM Amplitude 2 mm Heating range of ambient temperature to 83°C (181.4°F) Ramp-up time of <5 minutes from ambient temperature to 55°C (±3°C) Difference in the temperature detected by the internal sensor and the temperature of the sample liquid is approximately –2°C
Pipetting system	Syringe size 1 ml Pipetting range 5–900 µl
Capacity	Up to 12 samples per run. Up to 16 samples per run with iST-NHS and TMT16plex
Touchscreen	10.4" TFT Touchscreen, active area 211.2 x 158.4 mm, resolution 800*600 SVGA
USB stick	USB2.0 Compatible OS: Windows 8, Windows 7, Windows Vista, Windows XP (SP3 or later); Mac OS X 10.1 or later Operating temperature range: 0 to 35°C Operating humidity range: 10 to 90% (with no condensation) Storage temperature range: –20 to 60°C (–4 to 140°F) Storage humidity range: 10 to 90% (with no condensation) Formatting: FAT32
Handheld scanner	Imager Sensor: Wide VGA (752 x 480 pixels) Motion Tolerance: 25 IPS Light Source: Aiming: 650 nm VLD Print Contrast Ratio: 25% minimum Decode Capability: Reads standard 1D, PDF, 2D, Postal and OCR symbologies
Software	PreOmics protocols are preinstalled on the PreON. To learn more, please reach out to your local PreOmics contact. If unsure, get in touch at preomics.com/contact

10.5 Declaration of Conformity

Name and address of the legal manufacturer:

PreOmics GmbH
Am Klopferspitz 19
D-82152 Planegg/Martinsried

Phone: +49-89-2314163-0
Fax: +49-89-2314163-99
E-mail: info@preomics.com

An up-to-date declaration of conformity can be requested from PreOmics.

10.6 Waste Electrical and Electronic Equipment (WEEE)

This section provides information about disposal of waste electrical and electronic equipment by users.

The crossed-out wheeled bin symbol (see below) indicates that this product must not be disposed of with other waste; it must be taken to an approved treatment facility or to a designated collection point for recycling, according to local laws and regulations.

The separate collection and recycling of waste electronic equipment at the time of disposal helps to conserve natural resources and ensures that the product is recycled in a manner that protects human health and the environment.

Recycling can be provided by PreOmics upon request at additional cost. In the European Union, in accordance with the specific WEEE recycling requirements and where a replacement product is being supplied by PreOmics, free recycling of its WEEE-marked electronic equipment is provided.

To recycle electronic equipment, contact PreOmics Technical Support for the required return form. Once the form is submitted, you will be contacted by PreOmics either to request follow-up information for scheduling collection of the electronic waste or to provide you with an individual quote.

10.7 FCC Declaration

The "United States Federal Communications Commission" (USFCC) (in 47 CRF 15. 105) declared that the users of this product must be informed of the following facts and circumstances.

"This device complies with part 15 of the FCC:

Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.”

“This Class B digital apparatus complies with Canadian ICES-003.”

The following statement applies to the products covered in this manual, unless otherwise specified herein. The statement for other products will appear in the accompanying documentation.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules and meets all requirements of the Canadian Interference-Causing Equipment Standard ICES-003 for digital apparatus. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that the interference will not occur in an installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

PreOmics is not responsible for any radio television interference caused by unauthorized modifications of this equipment or the substitution or attachment of connection cables and equipment other than those specified by PreOmics. The correction of interference caused by such unauthorized modification, substitution or attachment will be the responsibility of the user.

10.8 Liability Clause

PreOmics shall be released from all obligations under its warranty in the event repairs or modifications are made by persons other than its own personnel, except in cases where the Company has given its written consent to perform such repairs or modifications.

All materials replaced under this warranty will be warranted only for the duration of the original warranty period, and in no case beyond the original expiration date of original warranty unless authorized in writing by an officer of the Company. Read-out devices, interfacing devices, and associated software will be warranted only for the period offered by the original manufacturer of these products. Representations and warranties made by any person, including representatives of PreOmics, which are inconsistent or in conflict with the conditions in this warranty shall not be binding upon the Company unless produced in writing and approved by an officer of PreOmics.

11 Appendix B – PreON Accessories

For more information and an up-to-date list of available protocols, visit <https://preomics.com/downloads>.

Product	Contents
PreON ⁸ Starter Pack, PreON	Instrument and 1-year warranty on parts and labor rotor adapters (12); spin column-based cartridge adapter (12); 200 µl Tips, racked (4x32)
Accessories	
iST PreON consumables	For 96 preps: 96 Disposable Rotor Adapters, 96 spin column-based cartridge adapters and 100 microcentrifuge tubes (1.5 ml); for use with the PreON
iST PreON Tips (1024)	Disposable Filter-Tips, racked; (8 x 4 x32)
Rotor Adapters (10 x 24)	For 240 preps: 240 Disposable Rotor Adapters and 240 microcentrifuge tubes (1.5 ml); for use with the PreON
Rotor Adapter Holder	Holder for 12 disposable rotor adapters; for use with the PreON
USB Stick	USB Stick; for use with the PreON
O-Ring Change Tool	O-Ring change tool for use with the PreON
O-Ring Set	Set of 10 O-Rings for use with the PreON

For up-to-date licensing information and product-specific disclaimers, see the respective PreOmics kit handbook or user manual. PreOmics kit handbooks and user manuals are available at www.PreOmics.com or can be requested from PreOmics Technical Services.

⁸ All configurations include: PreON robotic workstation for automated purification of peptides using iST Kits , bar code reader and 1-year warranty on parts and labor.

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