

vashe[®]
wound solution

Never Compromise

Safety or Efficacy





86%


WOCN reported on average **4+ departments** are involved in wound treatment¹



WOCNs reported an average of **4+ prescribers** are involved in wound treatment plan¹



Almost $\frac{1}{2}$ of WOCNs have **4+ cleansers** to choose from for wound treatment¹



Evolving Practice means Never Compromise

Wound bed preparation and cleansing technologies have significantly evolved in the last decade. Newer products now take into consideration the need for a cleanser that can disrupt microbial colonies while maintaining the integrity of key wound healing cells.

As a result, expert guidelines now include evidence-based recommendations that can lead to improved outcomes and cost savings without compromising safety and effectiveness.

Following evidence-based consensus guidelines



Limited evidence exists on the ability of saline to address high levels of microbes or bacteria present in chronic wounds, while other antimicrobial preservatives present in cleansers could compromise wound healing. Using a cleanser that can remove or disrupt microbial colonies is a critical component to wound healing.^{2,3}

The ideal cleansing solution should balance the need for removal of microbial colonies while avoiding damage to key wound healing cells. Cytotoxicity to these key cells; fibroblasts, keratinocytes, vascular and endothelial cells should be considered as they are required to heal the wound.²⁻⁵

Traditional cleansers such as hydrogen peroxide, traditional sodium hypochlorite (e.g. Dakin's solution), povidone-iodine and chlorhexidine are proven to be cytotoxic to necessary healing cells, which leads to “hard-to-heal” wounds burdening healthcare. A wound cleanser for frequent use should both disrupt and remove germs/microbes and preserve wound cells to promote a healthy wound environment.²⁻⁵

vashe[®]
wound solution

**Never
Compromise.**



Meets all new consensus guidelines requirements




Allows you to standardize your practice



Enables continuity of care that can lead to better outcomes.

Replicating the body's natural response to invading pathogens:

After a pathogen enters the body, neutrophils are quickly deployed from the bloodstream to respond. Through a complex biochemical pathway, *Pure* hypochlorous acid (pHA) is produced naturally by the human body to neutralize invading pathogens.



Proven efficacy

A recent prevalence study confirmed that almost 80% of chronic wounds contain high levels of microbes.^{6,7} Those microbes are thought to be the root cause of ~80% of all infections in humans and responsible for potentially delaying healing in 60% of chronic wounds.^{8,9}



Pure hypochlorous acid (pHA) has the ability to disrupt microbial colonies after short exposure¹⁰



No known clinically relevant resistance to **pHA**, non-mutagenic properties unlike other solutions (antimicrobials and antibiotics)



80-100 times more effective than sodium hypochlorite^{11,12}



Effective as a preservative against fungi, spores, viruses and multi-drug resistant bacteria

When a pHA-preserved wound cleanser is used in clinical studies, significant quantities of pathogens are mechanically removed from wounds, allowing the immune system to sustain the reductions.¹³⁻¹⁵

Organism	Time to kill	% Reduction
MRSA	15 seconds	99.999%
VRE	15 seconds	99.999%
<i>Escherichia coli</i>	15 seconds	99.999%
<i>Acinetobacter baumannii</i>	15 seconds	99.999%
<i>Bacteroides fragilis</i>	15 seconds	99.999%
<i>Candida albicans</i>	15 seconds	99.999%
<i>Enterobacter aerogenes</i>	15 seconds	99.999%
<i>Enterococcus faecium</i>	15 seconds	99.999%
<i>Haemophilus influenzae</i>	15 seconds	99.999%
<i>Klebsiella oxytoca</i>	15 seconds	99.999%
<i>Klebsiella pneumoniae</i>	15 seconds	99.999%

Organism	Time to kill	% Reduction
<i>Micrococcus luteus</i>	15 seconds	99.999%
<i>Proteus mirabilis</i>	15 seconds	99.999%
<i>Pseudomonas aeruginosa</i>	15 seconds	99.999%
<i>Serratia marcescens</i>	15 seconds	99.999%
<i>Staphylococcus epidermidis</i>	15 seconds	99.999%
<i>Staphylococcus haemolyticus</i>	15 seconds	99.999%
<i>Staphylococcus hominis</i>	15 seconds	99.999%
<i>Staphylococcus saprophyticus</i>	15 seconds	99.999%
<i>Streptococcus pyogenes</i>	15 seconds	99.999%
<i>Staphylococcus aureus</i>	15 seconds	99.995%
<i>C. difficile</i> endospores	15 seconds	99.93%

Proven safety



Based on years of clinical experience, evidence and extensive testing, Vashe helps to accomplish the goals of wound bed preparation and has proven to be:



Non-irritating



Safe for key cells



Safe around
mucous membranes



Has no known
contraindications



FDA cleared



Safe for key cells

Animal Model	Results
Eye Irritation (Rabbit)	No ocular irritation
Skin Sensitization (Guinea Pig)	No skin sensitization, no delayed-contact hypersensitivity
Primary Dermal Irritation (Rabbit)	No dermal irritation, no erythema or edema
Acute Oral Toxicity (Rat)	No oral toxicity (LD50>5g/kg)
Cell-Based Assay	
Bacterial Mutagenicity	Non-mutagenic
Cytotoxicity	Biocompatible with fibroblasts and keratinocytes

Wound Irrigant	Results	% Cell Survival (Fibroblasts & Keratinocytes)
Hypochlorous Acid (@ 4x the normal % of Vashe Wound Solution)	Pass	> 75%
Saline (0.9% NaCl, pH 5.0)	Pass	> 75%
Dakin's Solution (0.25%)	Fail	< 25%
Dakin's Solution (0.5%)	Fail	< 25%
Chlorhexidine gluconate (4%)	Fail	< 25%
Hydrogen peroxide (3%)	Fail	< 25%
Povidone iodine (7.5%)	Fail	< 25%
Povidone iodine (10%)	Fail	< 25%

Hypochlorous acid (at four times the normal percent of Vashe) is non-cytotoxic (grade 0), in contrast to other commonly used cleansers with significant cytotoxic effects¹⁸

A study was conducted in an outpatient wound center where Vashe was used for general cleansing on 31 patients. This study found that:

- **86% of chronic wounds healed**
- **Pain was reduced from 4.7 visual analog scale (VAS) to 0 at the end of the evaluation**
- **Odor was reduced from 4.58 VAS to 0 at the end of the evaluation¹⁷**

Addressing Patient Comfort




The importance of pH

Wound healing is optimal in slightly acidic environments where antimicrobial properties are higher. *Pure* hypochlorous acid has a pH between 3.5 and 5.5, which is favorable to wound healing environments that can aid in:¹⁸

- Optimal protease activity and oxygen release¹⁹
- Reduced toxicity of bacterial end products²⁰
- Epithelization and angiogenesis^{21,22}
- Increased macrophage and fibroblast activity^{21,23,24}

Many cleansing solutions contain toxic ingredients, such as sodium hypochlorite, and have a highly alkaline pH. An alkaline environment can allow pathogens to thrive and potentially impede the healing process.²⁵

Chlorine, Hypochlorous Acid, and Sodium Hypochlorite Abundance Based on pH + Relative pH of Wound Types and Various Solutions Used



vashe[®]
wound solution

A unique wound solution
manufactured at a pH of 5.5

Vashe in clinical practice

The effectiveness and versatility of a pH-A-based cleanser allows for streamlining and standardization. Vashe can be delivered through cleansing, packing, soaking, adjunctive debridement modalities (such as enzymatic and ultrasonic debridement), and in conjunction with Negative Pressure Wound Therapy with instillation and dwell (NPWT-id) for a variety of skin/wound conditions.



Venous Leg Ulcers



Pressure injuries
(Stage I-IV)



Burns



Post-surgical wounds



Diabetic Foot Ulcer



Pediatrics



Trauma

Unmatched Support



With so many potential variables, consistent wound cleansing seems impossible. When asked, almost **half** of clinicians said evidence-based practice is only attainable with continuous education and compelling evidence.¹



Standardize

With an **evidence-based custom protocol** that meets your facility's needs & latest consensus guidelines.



Reduce SKUs

With a versatile solution that's appropriate for use in all departments



Empower Staff

With education your way: webinars, in person, or via many peer-reviewed publications

The Right Partner
& Product



To discuss your custom protocol contact your URGO Account Manager or call 1-855-888-8273



Economic Benefits

A study was conducted among 24 patients with complex wounds including: necrotizing fasciitis, trauma, and pressure injuries. Negative Pressure wound therapy with instillation was initiated with Vashe vs. saline.²⁶

Conclusions:



7 Day reduction
in healing time

\$141,280

Cost savings per patient

Another study was conducted consisting of 17 adult patients with complex wounds, of multiple etiologies. Patients were randomly assigned to receive either Vashe or saline irrigation during low frequency ultrasonic debridement.²⁷

Conclusions:

55% Less complications

\$3K-\$33K Cost savings per complication

Vashe ordering information

Bottle Size/ Pack Size	Pour-Top Bottle Catalog Code	Instillation Bottle Catalog Code
4.0 fl. oz. (118 ml) Bottles/24-Pack	00312	Not available
8.5 fl. oz. (250 ml) Bottles/12-Pack	00313	00316
16.0 fl. oz. (475 ml) Bottles/12-Pack	00314	00317
34.0 fl. oz. (1 liter) Bottles/6-Pack	00322	00323



Vashe is also available over the counter! Visit vasheotc.com to learn more and order.

References: 1. Urgo data on file. 2. Murphy C, Atkin L, Swanson T, Tachi M, Tan YK, Vega de Ceniga M, Weir D, Wolcott R. International consensus document. Defying hard-to-heal wounds with an early antibiofilm intervention strategy: wound hygiene. *J Wound Care* 2020; 29(Suppl 3b):S1-S28. 3. European Pressure Ulcer Advisory Panel, National Pressure Injury Advisory Panel and Pan Pacific Pressure Injury Alliance. Prevention and Treatment of Pressure Ulcers/Injuries: Clinical Practice Guideline. The International Guideline. Emily Haesler (Ed). EPUAP/NPIAP/PPPIA: 2019. 4. Eriksson, E, Liu, PY, Schultz, GS, et al. Chronic wounds: Treatment consensus. *Wound Rep Reg*. 2022; 30(2): 156- 171. doi:10.1111/wrr.12994 5. International Wound Infection Institute (IWII) Wound Infection in Clinical Practice. *Wounds International*. 2022. 6. Malone M, Bjarnsholt T, McBain A.J., James G.A., Stoodley P., Leaper D., Tachi M., Schultz G., Swanson T., Wolcott R.D. The prevalence of biofilms in chronic wounds: A systematic review and meta-analysis of published data. *J. Wound Care*. 2017;26:20-25. doi: 10.12968/jowc.2017.26.120 7. Mendoza R.A., Hsieh J.-C., Galiano R.D. Wound Healing—Current Perspectives. IntechOpen; London, UK: 2019. The Impact of Biofilm Formation on Wound Healing. 8. Percival SL, Vuotto C, Donelli G, Lipsky BA. Biofilms and wounds: an identification algorithm and potential treatment options. *Adv Wound Care (New Rochelle)*. 2015;4(7):389-397. 9. James GA, Swogger E, Wolcott R, et al. Biofilms in chronic wounds. *Wound Repair Regen*. 2008;16(1):37-44. doi: 10.1111/j.1524-475X.2007.00321.x. 10. Robson MC. Treating chronic wounds with hypochlorous acid disrupts biofilm. *Today's Wound Clinic*. 2014;8(9). 11. Block SS. Disinfection, sterilization, and preservation. Philadelphia: Lea & Febiger; 2000. 12. Wang L, Bassiri M, Najafi R, et al. Hypochlorous acid as a potential wound care agent: Part I. Stabilized hypochlorous acid: A component of the inorganic armamentarium of innate immunity. *J Burns Wounds*. 2007;6:65-79. 13. Vashe Wound Solution data developed from USP 51 Antimicrobial testing. 14. Bohn GA, Champion S, Eldridge K. Can the use of hypochlorous acid change your dressing selection? Poster Presentation: Symposium for Advanced Wound Care: 2013; Orlando, FL. 15. Nerandzic MM, Rackalyte E, Jury LA, et al. Novel Strategies for Enhanced Removal of Persistent Bacillus anthracis Surrogates and Clostridium difficile Spores from Skin. *PLoS One*. 2013;8(7):e68706. doi:10.1371/journal.pone.0068706. 16. Data on file with Urgo Medical North America. 17. Niezgodna JA, Sordi PJ, Hermans MH. Evaluation of Vashe Wound Therapy in the clinical management of patients with chronic wounds. *Adv Skin Wound Care*. 2010;23(8):352-357. 18. Nagoba BS, Suryawanshi NM, Wadher B, et al. Acidic environment and wound healing: A review. *Wounds*. 2015;27(1):5-11. 19. Vermeulen H, van Hattem JM, Storm-Verstout MN, et al. Topical silver for treating infected wounds. *Cochrane Database Syst Rev*. 2007;24(1):CD005486. 20. Gethin G. The significance of surface pH in chronic wounds. *Wounds UK*. 2007;3(3):52-56. 21. Hunt TK, Hopf HW. Wound healing and wound infection: What surgeons and anesthesiologists can do. *Surg Clin North Am*. 1997;73(3):587-606. 22. Nagoba BS, Gandhi RC, Wadher BJ, et al. Microbiological, histopathological and clinical changes in chronic wounds after citric acid treatment. *J Med Microbiol*. 2008;57(5):681-682. 23. Thomas S. *Wound Management and Dressings*. London, UK: Pharmaceutical Press; 1990. 24. Molan PC. Re-introducing honey in the management of wounds and ulcers-theory and practice. *Ostomy Wound Manage*. 2002;48(11):28-40. 25. Ortega-Peña, Silvestre & Hidalgo-González, Christian & Robson, Martin & Krötzsch, Edgar. (2016). In vitro microbicidal, anti-biofilm and cytotoxic effects of different commercial antiseptics; Antibiofilm effects of different commercial antiseptics. *International Wound Journal*. 14. doi:10.1111/iwj.12625. 26. Gallagher KE, Alberto EC, Mallow PJ, Hermans MH, Cardenas L. A Retrospective Health Economic Analysis of a Stable Hypochlorous Acid Preserved Wound Cleanser Versus 0.9% Saline Solution as Instillation for Negative-Pressure Wound Therapy in Severe and Infected Wounds. *Cureus*. 2022 Apr; 20:14(4):e24321. doi: 10.7759/cureus.24321. PMID: 35607546; PMCID: PMC9123387. 27. Mallow PJ, Hiebert JM, Robson MC. Cost-Effectiveness of Hypochlorous Acid Preserved Wound Cleanser versus Saline Irrigation in Conjunction with Ultrasonic Debridement for Complex Wounds. *JHEOR*. 2021;8(2):76-81. doi:10.36469/001c.28429